

Product datasheet for **SC213159**

GRAMD2B (NM_001146319) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	GRAMD2B (NM_001146319) Human 3' UTR Clone
Vector:	pMirTarget (PS100062)
Symbol:	GRAMD2B
Synonyms:	GRAMD3; NS3TP2
ACCN:	NM_001146319
Insert Size:	1438 bp



[View online »](#)

Insert Sequence:	<p>>SC213159 3'UTR clone of NM_001146319 The sequence shown below is from the reference sequence of NM_001146319. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site</p> <pre>GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC TTACAGAAGTTGCTTGAGAATGGTGACTGATCGACCAGATTGCTTGGCCATCGGAATACCTCATGTTT CCCTTTGAAGAAGGTGCTTCTGAGGCGTTTTGTTTGTGAGTGCACCCTGCTGGTCAGAGGTGCAAGCAGA TGAGAATCCAGACATTGCATGTGGAGGTCTCCAGCTCAGGAAAGTGAGGAGGGAAATAATTTTGGTTCT TGTGCAATAAAAAGTTGACCTTGACTCTCTGAGGAAGATTTTGTGCTTTTGCCTGAAGAAAAACAGACCC ATCTCTGGAGGTCTCAGGAAGGGCCAGCGAACACTCTCTGGATAATTACCACGATGGCGTCAGCA AACACTCCACCCTGTGCCTTTTTAGTCCTTCCCGCCCTCTGCCTCTCCCTTACACCCTCTTAACGAC TTTCAAATAAAGGATACATCATATACTGACAACTCAATGTGGTCTTTCAAGAATTAGCCATGAGTC TCAAAAAGGCAATAAATGGCTCTAAGTGGACAGGTTTGTCTCAAACAAGTAACATCTACATTTTGTCTT TTTTTTTTTCAGTTCTCCTGTTATGTTCTGGTTGAAATCACCTGTGTGCTTAATTTCTCAATTCCTTTT TGGCAAGAATATCAAGCAAGGTGAATTTAACATTATGTTTATGTTTGTGTTTGTGCTGTAACATAATAG TTAATTGGACTGATTCTTACCCAGCCCTGGTCAAGAATCTGTGAGGCATGTGACTGAAGTACTAAATTA AACTTATTTTGAACCAACCTAATTTTAAAGCCAAAAGGTGTAATAGTGATTTAATACAGGATGAAAA ACACTGAATTTTTAAGACTGTAGGTGGACTATGTTAGTAGTTTTCAAGCAGGATGTCTGTATTACAGCAT TCAATAATGCTAAAATCCCTTTCAGCATGAAATTTGTATGTTTTATCCTTTGCTGACTAAAATAAAT AACTGGTGGTTTTGCTAACTATTTACCATATGGGTTGGCTGGTTTTTATTAATAATTGTTTAGGGTATCA TAAGATCTGTAAGTACAAAGGATATTCTGTTTCTATCTGCAATGCATTTGATAATCATTCTATGAAAT GTATTTTATTTAATCAGATAAGCTAATAAGCGAATTGGTTACATCGTTTGTATCTGTTGGCCTAGTGGA CAGCTGTGCCTGGTGCACCTCTGGCTTTTAATAAATACTATTGGTTGAAACACCCTGCTAACACTGTG TTTTGTGAGTTTCGTAATTCCTTTAACCAAGACTTAACATGCTTTTTAACTCTTTTTCCACCTGCCAT CTACTACTACCAACTCACTTATACTTCTTCTCCCAAGTAATCATTTCAAAACTAAAAAAGGATTTAA GTTCTTGGGTTACAACGGTTTAAAGGAAAGAAAATGAACTTTGAAGTCAAAAAAAAAA ACGCGTAAAGCGCCGCGCATCTAGATTCAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA CGAGATTTGATTCCACCGCCGCTTCTATGAAAGG</pre>
Restriction Sites:	Sgfl-MluI
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	NM_001146319.3
Locus ID:	65983
MW:	55.4