

Product datasheet for **SC213182**

PKIB (NM_032471) Human 3' UTR Clone

Product data:

Product Type: 3' UTR Clones
Product Name: PKIB (NM_032471) Human 3' UTR Clone
Vector: pMirTarget (PS100062)
Symbol: PKIB
Synonyms: PRKACN2
ACCN: NM_032471
Insert Size: 1208 bp
Insert Sequence: >SC213182 3'UTR clone of NM_032471
The sequence shown below is from the reference sequence of NM_032471. The complete sequence of this clone may contain minor differences, such as SNPs.
Blue=Stop Codon **Red**=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAACGATCGCC
TTGGAAAAGCCTCAAAATGAAGAAAATGAAGGCTCATAATCTATCAAGAGTGCTGAATTTCTGCATGT
TGAAAGACTTAGTGGTTCTGTTTTCTTGAGACATTTAATCTGGTGGTAACTGGTAACATTGCAGCCC
TAAGCAGCATGTGTATATTAGATAAATTGTGTTGTGATGCTACTCACTTTGATTGCAATGATGATGTCCA
AGGTAAGCTATTAAGGAGGTTACTTCCAAATCGCACTGAAGGAAAAGTTAAGAATAATACATGAT
CACAGAAATGCATACCACTGTCTGTAAACCAACAAAATCACTGTTCTTTTTGGATTTATTTAGCCT
GATGTATTTTTAATCAATTTTTATGGTATGGGCAAAATCATTCTTGGTAAATGTAATCAACATGAT
TGATTTAAACTTCATGGAATTTGTAGAAAATTATGGACATTTTTGGTGAGAAAGAACAATAGTCAAAA
CTCACATGGATAGAGTGTGTTTTGTTTGGCCAAAATGCCCGAGCTTTTCCCAAACCTCAAAAACG
TCTTGAAAAATGTAAAAGTTTGATAACAGAAACATCTTTAGGATATTTTTGTCTGACATATTTGCT
TCTAGTATGTGCCTACTGTGATTTTTTCATGTGGAAAATGCAAAATTTGTAACAAAATGGTTATATGG
AACATGCCTATTAATGAATTTACTATCTCCCTAACTTTGGTCTGTGTATGTGTGTGTTTTACTT
TAATATGAATTATACAAAATACTAGTTGTTTTACTCTCTTTTTCTTATTCTTAGGGCTTTTGTGTATG
TCTGACTGTTTTTAAATAAATCCTCAGCAATGCAGACCTTAATTTTTATTTTTTAAAGTAGCTA
ACATAGCAGTAGGCATTAAGCATTTAGTCAATGATATTGGTAGAAAATAGTAAAATACATCCTTTAAAT
ATATATCTAAGCATATATTTTAAAGGAGCAAAAATAAACCAAGTGTTAGTAAATTTTATTTTCA
GATATTTTAGAAAAATAAGAAATCTGAAGTTTTAAAAATGTCAGTAATTAATTTTATTTTCA
GAAATATATGCATGCAGTTATGTTTTATTTGATTGTTGACTTAGGCTATGTCTGTATACAGTAACCAAA
TAAACTTTTCACTATTAAGAGATTTCTACTGA
ACGCGTAAGCGGCGCGGCATCTAGATTGAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTTTCGATTCCACCGCCGCTTCTATGAAAGG
```

Restriction Sites: SgfI-MluI



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	NM_032471.6
Summary:	This gene encodes a member of the cAMP-dependent protein kinase inhibitor family. The encoded protein may play a role in the protein kinase A (PKA) pathway by interacting with the catalytic subunit of PKA, and overexpression of this gene may play a role in prostate cancer. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene. [provided by RefSeq, Jul 2012]
Locus ID:	5570
MW:	47.2