

Product datasheet for **SC213190**

HIP1R (NM_003959) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	HIP1R (NM_003959) Human 3' UTR Clone
Symbol:	HIP1R
Synonyms:	HIP3; HIP12; ILWEQ
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_003959
Insert Size:	1238 bp
Insert Sequence:	>SC213190 3'UTR clone of NM_003959 The sequence shown below is from the reference sequence of NM_003959. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAACGATCGCC
ATCTACCCAGCTCAACTCGTGAAGTACTAGGCCCCCAGGGGTCCAGCAGGGTGGCTGGTGACAGGCCCT
GGGCTCTGCAACTGCCCTGACAGGACCGAGAGGCCTTGCCCTCCACCTGGTGCCCAAGCCTCCCGCC
CCACCGTCTGGATCAATGTCCTCAAGGCCCTGGCCCTTACTGAGCCTGCAGGGTCTGGGCCATGTGG
GTGGTGCTTCTGGATGTGAGTCTTATTTATCTGCAGAAAGAACTTTGGGGTGCAGCCAGGACCCGGT
AGGCCTGAGCCTCAACTCTCAGAAAATAGTGTTTTTAATATTCTCTTCAGAAAATAGTGTTTTTAAT
ATTCGAGCTAGAGCTCTTCTTCTACGTTTGTAGTCAGCACACTGGGAAACCGGCCAGCGTGGGGCT
CCCTGCCTTCTGGACTCCTGAAGGTCTGGATGGATGGAAGGCACACAGCCCGTGCCGGCTGATGGGAC
GAGGGTCAGGCATCCTGTCTGTGGCCTTCTGGGGCACCATTCTACCAGGCCCTCCAGCTGCGTGGTCT
CCGCAGACCAGGCTCTGTGTGGCTAGAGGAATGTGCGCCATTACTCCTCAGGCCTGGCCCTCGGGCCT
CCGTGATGGGAGCCCCCAGGAGGGGTGAGATGCTGGAAGGGGCCGCTTTCTGGGAGTGAGGTGAGAC
ATAGCGGCCCCGGCGCTGCCTTCACTCCTGGAGTGAGTTTCCATTTCAGCTGGAATCTGCAGCCACCC
CCATTTCTGTTTTCCATTCCCCGTTCTGGCCGCGCCCACTGCCACCTGAAGGGGTGGTTTCCAGC
CTCCGGAGAGTGGGCTTGGCCCTAGGCCCTCCAGCTCAGCCAGAAAAAGCCAGAAACCCAGGTGCTG
GACCAGGGCCCTCAGGGAGGGGACCCTGCGGCTAGAGTGGGCTAGGCCCTGGCTTTGCCGTCAGATTT
GAACGAATGTGTGTCCTTGAGCCCAAGGAGAGCGGCAGGAGGGGTGGGACCAGGCTGGGAGGACAGAG
CCAGCAGCTGCCATGCCCTCCTGCTCCCCCAGCCCTAGCCCTTTAGCCTTTACCCCTGTGCTC
TGGAAAGGCTACCAAATACTGGCCAAGGTGAGGAGGAGCAAAAATGAGCCAGCACCAGCGCTTGGCTT
TGTGTTAGCATTTCTCCTGAAGTGTCTGTTGGCAATAAAATGCACTTTGACTGTTTGTGTCA
ACGCGTAAGCGGCCGCGCATCTAGATTGCAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG

```



[View online »](#)

Restriction Sites:	Sgfl-Mlul
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_003959.3</u>
Summary:	Component of clathrin-coated pits and vesicles, that may link the endocytic machinery to the actin cytoskeleton. Binds 3-phosphoinositides (via ENTH domain). May act through the ENTH domain to promote cell survival by stabilizing receptor tyrosine kinases following ligand-induced endocytosis.[UniProtKB/Swiss-Prot Function]
Locus ID:	9026
MW:	44.2