

Product datasheet for SC213199

PHLDB1 (NM_015157) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	PHLDB1 (NM_015157) Human 3' UTR Clone
Symbol:	PHLDB1
Synonyms:	LL5A
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_015157
Insert Size:	1235 bp

Insert Sequence: >SC213199 3'UTR clone of NM_015157
 The sequence shown below is from the reference sequence of NM_015157. The complete sequence of this clone may contain minor differences, such as SNPs.
 Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAACGATCGCC
GCTGAGGGCTACACTCAGTTCATGAACAACTGCCGTGGGCTCCTGGCAGAGCACAACTGGGGCTTTT
GTATAAGAAGACTTTAATATTCTGTAAGGAGCTTGGTCTGTGAGTTTCTGGGCTCTGGCCTCCTGAAG
AACCAGCCAGAAGAAGAAAAGTAGAGGTGGCTTTGCTGCCTCCTGGGAGCCAGAACTTGCAGTAACCC
TTTAGGGTCTGCCCCAGGCCAGCCAGGGCTGAGGAGCTGTACAGAGAGGGCCTCAGCTCTGACCTG
ACACCTGCTCTCCCAGCCTGTTTTCTCTTTCTAAAAGACAAATTATGGTACCATAAGCTGCCAAAGA
TCCCTCTGCTCAGACCCCTTTGCCAGGGGCTTTGGGGCTGAGCAGAGCCACATCCAGAGTGGGGT
AATAGCTCAGGCGGCCGCTTCCCATTTCTCAAACCCGCTCTGCCCCATTGTTCTCCTTTCCCTTATA
CTTTTTATTACCTTGCTCAAGGGCCAGAGATCTCAAGTGTCAACCTTGAGGTCCCAGCTCCATCCCTA
GTTGCAGACTCATCACCATTGTTACCATAGTACTGCTTCATTGCCATGGTTACATACTAATTGTGCA
GCTCTGTGGCCAGCCCACTGCTTCAGCTGTGGGCCATCTGAGGTACGTGCCATCATCTCTCCAGCCC
AGGCCCCCTGGGCATCTCATGCTGGGGGAAGGGACTGAATACCTTTTCCCTTCCCCTGCCTGTGTCTT
CAGCCCTGATGCACAGGCTGCCAGCCCCCAGTCCAGCCCTCTCCCTTCCACTGGTGCCTTGCTTAGAG
CCAGAAGGGATGAAGCCGGGGATCTATGGAACAGAGGAGGAGCGATGCAGTTGGGAGAGGAAGCTAGA
AGGGTTATGGTTGGAGTTCTGTACAGTGTGAGTTTCCGACAGGGAAGAGGATTCTCCAATGCTCCT
AGAGAGAAAGCCTGAGCAGGAGATGATGCAGCAGAGGGGAAGGGCCCTGTGGTGGCGCCGCCCTTCTT
CAGCCTCCGAAGGGTGTGGAATGGAGAGTGGAGGACCAGGCCTCCAGCTGTCTGGCCTCGCCCTTCA
CGCCTTAACACTAAGCCACCTCCCCTGCTCTCCTTCCAGCATTGAGCCCTTGTTGCCTGGGCCAG
GCTGGGGGTTTTAGTATTTGTAAGCATTTAGCAGACAATAAAGCCTTTGGACTACGGAA
ACGCGTAAGCGCCGCGCATCTAGATTGGAAGAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
  
```


[View online »](#)

Restriction Sites:	Sgfl-Mlul
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_015157.4</u>
Locus ID:	23187
MW:	44.7