

Product datasheet for SC213205

ISM1 (NM_080826) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	ISM1 (NM_080826) Human 3' UTR Clone
Symbol:	ISM1
Synonyms:	bA149I18.1; C20orf82; dj1077I2.1; ISM; Isthmin
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_080826
Insert Size:	1222 bp
Insert Sequence:	>SC213205 3'UTR clone of NM_080826 The sequence shown below is from the reference sequence of NM_080826. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAACGATCGCC
AAGCAGTTCCAAGAGGCCAGGGAATATTAAGAGACTGGGATGAGGTGGAGGACGCTGCCTCTGTTCT
GGAGCACACACGTGCTGCACTGACGTGCCGACTGGCGCCGAGACCTTCATAGCTGCGGTCTGTATATT
TGTATATACCACATGAGTATTTCTCATACATTACGCTAGGGCGTGTGCCACGCCAGGGGACTGCCTT
GTGAAGCCGCCCTCGCCATCTGCAGAGCTCCTTGAAGTCCCCCTGGGGAGCGATGTGGCAGAAGGAT
GGGGACAACCTTGAAGCCAGAAGAAGAACTGGAAGCCACAGTGGGTGCGACTCAATTCACACCCGAT
CCAGAGTTTCAAAGAGAGGCAAAGGGGAAAGAGACTGAGGTTGTAAACGTTATAAGCAGTTTTATAT
ATAACTTATTTAATAACAATGTGACTTAATTAAGCGTAACCTTTCTCTGGAGTTGTGGTGAAACTAAT
CACGTCTGTGAGAGATCAGAAAGAAAGAGACTTAGGGAAGTGGAAGAGAAAGGGAATTTTGAATTTAT
TTCTTTAAAAATAATGCAATGGAATATATCAAAACATGTAAACGCCACCTTAAACCAATGTTATTT
GGTCATGAGGCACCTTGCTGGAGTCTCAGATTCCAAAAGTCTCTTCTTACAGCTGGGTAGGAATATGA
TATTTTAGGGACAAAGCTGAGGACTGGTTTTAAATAGGCTTTAAATAAAAGATCAATATTATCATAAT
GCTATCATTCTGCTAAACGGCCCCAAAACAGTAGAATTTCTGCTCATGTCCTAGCAGGTTCAAGAGACT
GCAGCCAAGTTCAGATGTAAAAACAAGAGTAGCACTTTTCCAAAGGAAAACAACAAAACAAATGGGAA
AAAGATAATGGACCGCATTTACCTATTTTAATACTATTTTAACAATTTTTTCATCTACCAATATATCC
CCAATAAATAAATAAAAGGGGGGAGGGTCAATCTGGGGAATCTTAGTTTTTTATGTTTTAAGAAAA
CAAAAAAACTGCATTATTTTTGTAAAGTATTTATTGAGTCACGGATTATTGTGCATCAAGCAATTGTT
AATATGACCTGGTCTATGGGTAGAACTTAGGAAAAATAAAGTTGGTTCTTATTCAATATTTTACTTT
GCAAAATTCTAGTAAAGAGAGTATATAATAAAATCATAATAAAAGGTG
ACGCGTAAGCGGCCGCGCATCTAGATTCAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
  
```


[View online »](#)

Restriction Sites:	Sgfl-Mlul
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_080826.2</u>
Summary:	Acts as an angiogenesis inhibitor.[UniProtKB/Swiss-Prot Function]
Locus ID:	140862
MW:	47.6