

Product datasheet for SC213206

ANKRD13B (NM_152345) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	ANKRD13B (NM_152345) Human 3' UTR Clone
Symbol:	ANKRD13B
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_152345
Insert Size:	1209 bp
Insert Sequence:	>SC213206 3'UTR clone of NM_152345 The sequence shown below is from the reference sequence of NM_152345. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAACGATCGCC
ATCCTGAGGCTCTCACTGACCGAGCAGTAGCGCCCCCTGCCGGGACCCCTCGCCAGCGCCACGCGGCCA
CGCCAGGGCCAGGAGCCAGACAAACCCGGCCTGCGCGCCTGCAGAGCGCGGCTGGAGACTGGAGCC
ACCGCCTCGCGGGTGCAGCAGCACAGCAGGCACGGTTCGGGGAGGGATTTCGGCATGGCCGCGGGTACC
TTCCCAGGCCAGGGCCCTGGAGGCAACTGGCAGGCCTGGTCCCCCTGCTTTGCTGTATTCTGATTCCC
CAACCCGCTCCCCTGGGCTCAGATCTGTCTGTCTAGGGCGGAGCCAGCGGTCTGAGGGGAGATG
AATCCTTAGAGGAGCGCTGTCCCTATCCCTTGCTCCTTGCGCCGGTCCGCATTCTGTTCACTAACCC
CACTCCAGGACTGCCCCTGAAGCCTTTGCATCTCTGCTCTTCACTCCTGGGGGCGAGTGAGTCCC
CGCGTGCTGTGCTACTGCTTTGTTCCAGACACAAACCAGCAGTCTAGGGCCAGCCCCCCCCACCC
CGGCATTCAGCGTCAAGTGCACTTAGCGGGGTACCCGGCTCCCCAGCCCCACACCCGCTCCTGGGTC
TCAGGGTGGTTCAGCTTCTCCTTGGGCGCCAGAAGTTGGAGTCCCATCCCCAAGGCACATTTTTTT
TTTTTTTTTTTTGTGATCAGGGAGTGTGTCACGTAGCCCCCTGGCCTGTGCAAGCCCTGACTCCCTC
ATGGTGCCTCGGAGAGTGGGAGCATATTGGGCTGGGGTAAGCACTAGACCCAAGTAGACTGGACACAA
AGGGCTCGCCAGGGCCCTGGCGCCACCCCAACCCCTTCCACCAAGCTGCTGCTAGCCTCTGTGGTTGT
ACATCCCCTTGCCCCACACGGAGACTGACTCTAAAACCTTCAATCAATGGTGCTAACCCCCGGCTC
TCCCCTGCCCCACCTCACCCACCCAGAGAAGCACAGACCCCGCCAGGGGCGAGGGGCCACCGCACACCC
TTGTCCCGGGCCTGTCTGGGACTGGCCTTCCCGGCTCAGCCAGTGAGGCTCAGAAAGGACACAAAGAGG
GATGGAAGAAAAGACAAAGAGAACTGTTCTCCCAACCCCTTCCCTGATGCCAGGGGCACCAAGACTG
ATTCTGAGGCACAAATAAAGAGGCTTCATACCGGA
ACGCGTAAGCGGCCGCGCATCTAGATTCAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG

```

Restriction Sites: SgfI-MluI



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_152345.5</u>
Summary:	Ubiquitin-binding protein that specifically recognizes and binds 'Lys-63'-linked ubiquitin. Does not bind 'Lys-48'-linked ubiquitin. Positively regulates the internalization of ligand-activated EGFR by binding to the Ub moiety of ubiquitinated EGFR at the cell membrane. [UniProtKB/Swiss-Prot Function]
Locus ID:	124930
MW:	43.8