

## Product datasheet for **SC213209**

### NECAP2 (NM\_001145278) Human 3' UTR Clone

#### Product data:

**Product Type:** 3' UTR Clones  
**Product Name:** NECAP2 (NM\_001145278) Human 3' UTR Clone  
**Vector:** pMirTarget (PS100062)  
**Symbol:** NECAP2  
**ACCN:** NM\_001145278  
**Insert Size:** 1218 bp  
**Insert Sequence:** >SC213209 3'UTR clone of NM\_001145278  
The sequence shown below is from the reference sequence of NM\_001145278. The complete sequence of this clone may contain minor differences, such as SNPs.  
**Blue**=Stop Codon **Red**=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
CAGCCAGGCACAGGCTGGGTCCAGTTCAGCCTGAGCACGGTTTTCTCATGTGACTTCTGGGAAGGC
GCTCCCTCATCTGGGCCAAAGGAAGGAGGACGAAGCCCTCCTCAGCTGGCCTGTGTTGGGGCATGAAT
CTCTCCTCTCCTCTGTCTGGCTCTGTTGACAAACCGGCATGTTTGGCAGTAAATTGGCACCGTGTC
ACACTGTTTCTGGGATTCAAGTATGCAACCAGAACACAGGAGAAGAAAAGCTCCAGGATCCCTGTCCC
CATCTGCTCTTGTATGTGAGAGAGACTCTGAGACTTCTTCCATCGCAATGACCTGTATTAACACAAG
CCCCCAAGCAAAAAGAAGAGGTTGAGTTTGTGCCAGGATTAGATCAGCCCTTCCCAGGGTCTGCAGG
TGTACATGATCACAGTTCAGCGGGAGGCTTCCGTACCCACACTGGCTGTAGCCACTTCAGTCCATCT
GCCCTCCAGAGGAGGGGTTTCTCCTGATTTTTAGCAGGTTTAGAGGCTGCAGCTTGAGCTACAATCAG
GAGGAAAATTGGAAGGATTAGCAGCTTTTAAAAATGTTTAAATATTTTGCCTTGTCTAATGTGCTGATCC
GCACTAACTCATCTTTGCAAAAAGAACTGCTCCCTCGGCGTGCCCAAGCTGGGGCCTCTGAAGGGATTC
CTCACTGTGGGCAGCTGCCCTGAGCTTTCAGGCAGCAGTGTATCTCTGGCCAGTTGTCTGGTTTTCCAT
GTATTCTAGGCCAGGTAGGCAACACAGAGCCAAGGCGGGTGTGGAAGCCAGACGGAACAGTGTGGGG
CAGGAAGGTGGATGCTGTGTGTCATGGAGCTGTGGGAGTTGGCACTCTGTCTGCTGGTGGCCCTCTCGGC
TCACATGTTACAGTGCAGCTCCTGGCAGACTTGGGTTTTCTCTTTGGTGGTTTTCAAAGTGCCTTATC
TGCAAACAACCTTTTTCTCCTCAGGAACTGTGAATGGCTAGAAGAAGGAGCTCAGTAACTAGAAGT
CCAGGGTTGCTTGGTTTACTGGTTTATAAGAAATCTGAAAGCACCTCTGACATTCCTTTTATTAECTCA
CCTCTCAGTTGAAAGATTTCTTCTTTGAAAGGTCAAGACCGTGAAGTGAAGGAGTGTGGCCTTTTTG
CGGACCAGATTTTTAAGATAAAATAAATATTTTTACTTCTGTCA
ACGCGTAAAGCGCCGCGCATCTAGATTGAAAGAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
```

**Restriction Sites:** SgfI-MluI



[View online »](#)

<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
<b>Components:</b>	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
<b>RefSeq:</b>	<a href="#">NM_001145278.2</a>
<b>Summary:</b>	This gene likely encodes a member of the adaptin-ear-binding coat-associated protein family. Studies of a similar protein in rat suggest a role in clathrin-mediated endocytosis. Alternatively spliced transcript variants have been described. [provided by RefSeq, Feb 2009]
<b>Locus ID:</b>	55707
<b>MW:</b>	44.7