

## Product datasheet for **SC213236**

### Elastin (ELN) (NM\_001081752) Human 3' UTR Clone

#### Product data:

**Product Type:** 3' UTR Clones  
**Product Name:** Elastin (ELN) (NM\_001081752) Human 3' UTR Clone  
**Vector:** pMirTarget (PS100062)  
**Symbol:** ELN  
**Synonyms:** ADCL1; SVAS; WBS; WS  
**ACCN:** NM\_001081752  
**Insert Size:** 1237 bp  
**Insert Sequence:** >SC213236 3'UTR clone of NM\_001081752  
The sequence shown below is from the reference sequence of NM\_001081752. The complete sequence of this clone may contain minor differences, such as SNPs.  
**Blue**=Stop Codon **Red**=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAACGCATCGCC
GGGAAAGCTTGTGGCCGGAAGAGAAAATGAGCTTCCTAGGACCCCTGACTCAGACCTCATCAACGTTG
GTGCTACTGCTTGGTGGAGAATGTAACCCCTTTGTAACCCCATCCCATGCCCTCCGACTCCCCACCCC
AGGAGGGAAACGGGCAGGCCGGCGGCCTTGAGATCCACAGGGCAAGGAAACAAGAGGGGAGCGGCCAA
GTGCCCGGACCAGGAGGCCCCCTACTTCAGAGGCAAGGGCCATGTGGTCTTGCCCCCACCATCCC
TTCCACCTAGGAGCTCCCCCTCCACACAGCTCCATCTCCAGGGAACTTGGTGTACACGCTGGTGC
TCTTATCTTCTGGGGGAGGGAGGAGGGAAGGGTGGCCCTCGGGGAACCCCTACCTGGGGCTCCTC
TAAAGATGGTGCAGACACTTCTGGGCAGTCCCAGCTCCCCCTGCCACCAGGACCCACCGTTGGCTGC
CATCCAGTTGGTACCAAGCACCTGAAGCCTCAAAGCTGGATTGCTCTAGCATCCCTCCTCCTGGG
TCCACTTGGCCGTCTCTCCCCACCGATCGCTGTTCCCCACATCTGGGGCGCTTTTGGGTTGAAAAACC
ACCCACACTGGGAATAGCCACCTTGCCCTTGTAGAATCCATCCGCCATCCGTCCATTCATCCATCGG
TCCGTCCATCCATGTCCCAAGTTGACCGCCCGGCACCACTAGCTGGCTGGGTGCACCCACCATCAACCT
GGTTGACCTGTATGGCCGCTGTGCCCTGCCCTCCACCCCATCTACACTCCCCAGGGCGTGCGGGG
CTGTGCAGACTGGGGTGCAGGCATCTCCTCCCCACCCGGGGTGTCCCCACATGCAGTACTGTATACCC
CCCATCCCTCCCTCGGTCCACTGAACCTCAGAGCAGTTCCCATTCCTGCCCGCCCATCTTTTTGTGTC
TCGCTGTGATAGATCAATAAATATTTTATTTTGTCTGGATATTTGGGGATTATTTTGTATTGTTGA
TATTCTCTTTGGTTTTATTGTTGTGGTTCATTGAAAAAAGATAATTTTTTTCTGATCCGGGGA
GCTGTATCCCCAGTAGAAAAACATTTTAACTACTCTAATAAATACTCTGGATGAAACACACCTTTTTTT
TTAATAAGAAAAAGAGAATTAAGTCTTCAGAAATGACTAATAAATGAAAAACCTTTAAAGGAAA
ACGCGTAAGCGGCCGCGGCATCTAGATTGGAAGAAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTTTCGATTCCACCGCCGCTTCTATGAAAGG
```

**Restriction Sites:** SgfI-MluI



[View online »](#)

<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
<b>Components:</b>	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
<b>RefSeq:</b>	<a href="#">NM_001081752.3</a>
<b>Summary:</b>	This gene encodes a protein that is one of the two components of elastic fibers. Elastic fibers comprise part of the extracellular matrix and confer elasticity to organs and tissues including the heart, skin, lungs, ligaments, and blood vessels. The encoded protein is rich in hydrophobic amino acids such as glycine and proline, which form mobile hydrophobic regions bounded by crosslinks between lysine residues. Degradation products of the encoded protein, known as elastin-derived peptides or elastokines, bind the elastin receptor complex and other receptors and stimulate migration and proliferation of monocytes and skin fibroblasts. Elastokines can also contribute to cancer progression. Deletions and mutations in this gene are associated with supravalvular aortic stenosis (SVAS) and autosomal dominant cutis laxa. [provided by RefSeq, Aug 2017]
<b>Locus ID:</b>	2006
<b>MW:</b>	45.1