

Product datasheet for **SC213242**

FUSIP1 (SRSF10) (NM_006625) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	FUSIP1 (SRSF10) (NM_006625) Human 3' UTR Clone
Symbol:	FUSIP1
Synonyms:	FUSIP1; FUSIP2; NSSR; PPP1R149; SFRS13; SFRS13A; SRp38; SRp40; TASR; TASR1; TASR2
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_006625
Insert Size:	2000 bp

[View online »](#)

Insert Sequence: >SC213242 3'UTR clone of NM_006625
The sequence shown below is from the reference sequence of NM_006625. The complete sequence of this clone may contain minor differences, such as SNPs.
Blue=Stop Codon **Red**=Cloning site

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
TCTGCTTACTACACTTCAAGAAAGATCGAAAGCGGAAAAAGAACCAGAAAGGGCAGTTCAAGCGACC
AAAGGGTGGGTGGAAGGTGCTGCAATATGAATACTGTACGAATATTTTGAAGTAAAGCATTATAGGAC
AAATAAAGGAAATTCAGTTTGTACTGTGGAAGTAAATCCCTAAATATGAATAGGTTTATATTGATT
CATGGGTAAACAGGTCCATAATAAATTATTGAAAGTAAAGTAAATCAAGGAAGACAGCCATAG
TCTCTTACAGTGCCTCTGTTGGTCTGTCTCAAACTGAATTGGGTGGGAAAGGTATGGTCCAATATAAA
AGTTCCATTTTGGCATTATTGGCAAATCTGCCTTTGTTTATTTTGGTGCCAGTGTCTGCTTAAT
CATTTGCTTTGTTGGCATCTGTGTTTATTACTGTACACCACATGCAGTTTACATCTGTCTTAAGTAC
TCCTTCCAGGTAATTCGAATTATATTTGACATCCAGCTAAGAGGGCCCATCTCTCTACCTCTTTC
CTAGTCAGTATATTTCAGCAATATTTATTGAGCCCTTACTGTGGGCAAAATCATTGTACTGGATAATTGA
GAAAAATAGATAATTCCTTATTCAGTAAATGTCTACTGAGCACAATCTAGTGAATCATTACAGTATGG
CCTCATTTGTTTGTGAGGTGTGTTATTCATAACAATATTTTACACCATTCGTATCAATGTAATTATA
GAACACAATATACGATCAAGGATAAGTAATTGTGTGGTTATCTGCCATTTAAAGTATCCAGTATTTGA
TCACATTATTATAAATGAAGAAATGATTTAATCTGTAATAAACTGGTTTATTGTGCAGTGACTGTA
ATATACTAGAGTTATAAATGATTTTACTCTGCCTCACCAACACATGCTAGGATATAACCCCCAAA
TAAGTATTTAACTTTGCATTAGGTATAAAGGAGACTGGGTGCTATAATTAGATTATTTGAGGCAGACA
GAGAGCTGTTATCCTAACTGATTTAGTATGTTCTGTAATTGAGAAAAATGTTACCAAAATTATACTTTT
AGTGATTTACATGTACATTTTATAGGGGACATGTTCTGTGTATAGCGAATAAATACTTTTATAGTATC
ACAAAATGTTTGGATTCTTAGTTTCTTATTAGGATATGATGTTTCTTACCTATACATTGATTTTAT
CACATTGATTTTGTCTCTTTTGGCCACATACTTAAATTTTCCATATAAAATTAACAAAAATCT
GGAAGTGAATTTAAGCTTTTAAATTTGTAGCAATTTTGAAGTATAGCAGAGGTATTTTCAAGGGAG
GTGTGAAATGGAGGTATTATTGTTGATTTTATTTAGAAATTTGTTTGAAGTGGCCCAAGATCACA
CAGCAAGTGTGAGTCTAGGTTACAGCCCATATTTTCACTTGTTCATCTTTTGTACCATGTTTG
TTGTTTAAATCCAATGGGGAAGTGTGTTTATAGCTTTATTTTATGAGATAGTAACATTGATATAATGG
GATCTAGAATTTACTTAGAGTCAAAGTGAAGTATTAATTACAAAAAGTATTTTGTGTATACCGTTA
TACACTCATTAAACATAGGGGTTCTATGGACAGAGGTTCTATGATTTTTAAAGAAATATGTAGAGGCT
ATAGACTTTTTTGGCCTTACCTCCCTTCTCCATTTCTTTTGAAGTAAAGACATTTAAGGGGGAGAAT
TAAGGTGCATATTGAGTTGGGGTTTGTGTTAAAAAAATAGGTAATACATATATAACTACAATAACA
ACAGTTGGTAGTTCTCTAAATATTTAATTATGGTAAATCAGAAAAGGCCCTTTAATTTTGCATCTG
ACGCGTAAGCGGCCGCGGCATCTAGATTCAAGAAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
  
```

Restriction Sites: SgfI-MluI

OTI Disclaimer: Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).

Components: The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.

RefSeq: [NM_006625.6](#)

Summary:	This gene product is a member of the serine-arginine (SR) family of proteins, which are involved in constitutive and regulated RNA splicing. Members of this family are characterized by N-terminal RNP1 and RNP2 motifs, which are required for binding to RNA, and multiple C-terminal SR/RS repeats, which are important in mediating association with other cellular proteins. This protein interacts with the oncoprotein TLS, and abrogates the influence of TLS on adenovirus E1A pre-mRNA splicing. This gene has pseudogenes on chromosomes 4, 9, 14, 18, and 20. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Jul 2014]
Locus ID:	10772
MW:	77.3