

Product datasheet for **SC213243**

SLC25A45 (NM_001077241) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	SLC25A45 (NM_001077241) Human 3' UTR Clone
Symbol:	SLC25A45
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_001077241
Insert Size:	1245 bp
Insert Sequence:	>SC213243 3'UTR clone of NM_001077241 The sequence shown below is from the reference sequence of NM_001077241. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAACGATCGCC
TACGAATATCTCCTCCGCTGGTGGGGATGAGCCCTGCGGCAATGCCAGCAGCTCCCCATCAGGCCACG
GCCTGGAGGCCAGTTTGAGATTGGAGGCCAGTTGAAAGCTTGCAATCAGTGCAAGAGGCTCAGCCCT
TCCTAACCAAGGTGCCTCCCACCCGCGCAGATCTGGGCTGGGCAGACACCTGTGGGAGCCGGAAGCCAG
GGGGCCTGTGCAGCCTCCCTGTGTAGCTGGCCTTGACTCCTTTGCCTCCCACATCTGTGAACAGGGAG
CATGAGGCACAAGTGAGCTGGCAAGTGGTGTGTGACATCCCAGCTCCTGTCTGTGCCTTCACCTCT
TTTTTTTTTTTTTTTTTTTTTTTGTAGATGGAGTCTTTCTGTGTACCCAGGCTGGAATGCAGTGGCGTGA
TCTCGGCTCACTGCAACCTCCGCTTTCCACCTTCCGGGTTCAAGCAACTCTCCTGCCTCAGCCTCCCA
AGTAGCTGGGACTATAGGCGCCCGCCACCACCCAGCTAAATTTTGTATTTTAGTAGAGACAGGGT
TTCACCATGTGAGCCAGGATGGTCTTGATCTCCTGACTTTGTGATCCGCGCGCTCGGCCTCTCAAAGT
GCTGGGATTACAGGCGTGAGCCACCGCGCCGGCTGCCTTCACCTCTTAAGGAGCTCTGAGACTCCACT
TCTGAGAGTCCCTGCGGCTCCACCTCCCTGCCTTTCAAAGCTCTCTCCCCATGACCCAGGATAACC
CTATGTCTCCTCCCCAGAATCCTTCAGTGGCTCTCATCACCTTCAGGAAAAGCCCAAACCTCTCCAC
CCTCCAGGTTCTCCCTCCCACGAGGCTTTGTTCTCCTGGGTTGCTTCTTGACCTGAACAAGTTG
TGCTCTCATTGGCCGGGCTGGCCAGCAGTGCACAGTGCCTGGCAGGTTGACTCTACCATCCCCGGGGC
TGGCCCCGCTCTCCTCCGAGACCCAGGCTGAGCCAGTCCCCTCACCTTCTCTGACTTACCTCCCCAC
CTGAGGCTGACTTTGGGTTCCAGACACCTACCCACACACATGCCTTGATCATAGCACTTGCCTGCG
CTTCTTCAGAGTCATTAATTTGCTTCTCGGCTTCCCCACTGGAGTGTGAGCTGCCTGAGGTCAGGGATT
GCGCTTTGGATGGTTTCCAGCCTGAGCCTGGTGTGTAACAGACGTGTGAATAAATGCTCGTTAAATG
ATG
ACGCGTAAGCGGCCGCGCATCTAGATTCGAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
  
```


[View online »](#)

Restriction Sites:	Sgfl-Mlul
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_001077241.3</u>
Summary:	SLC25A45 belongs to the SLC25 family of mitochondrial carrier proteins (Haitina et al., 2006 [PubMed 16949250]).[supplied by OMIM, Mar 2008]
Locus ID:	283130
MW:	45.9