

Product datasheet for SC213248

SLC25A17 (NM_006358) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	SLC25A17 (NM_006358) Human 3' UTR Clone
Symbol:	SLC25A17
Synonyms:	PMP34
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_006358
Insert Size:	1234 bp
Insert Sequence:	>SC213248 3'UTR clone of NM_006358 The sequence shown below is from the reference sequence of NM_006358. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAACGATCGCC
ATGGGGCTGAAGCGTGACACCAACACAGACGCCTTCCCATGAAAAATCCGAAGATGCTCAAGAGG
GAGGTTTCTCTGAGTGAAGAGAAGTGATTCTCCCTTGACTCTGGCTCCTGCCACCACAAATGTTACC
CTCATTGGCTTGAAGCATCAAGGGTGACAGGGAGTATGGCCAACTGGACCTGTTGTCACCTTAAT
TGTCATGCTGGCTGGTTGGATTTGGGGTGGCAGTTGGACTAATGTGAAAAAACATTGCTGAAACCT
AAAAATGAAAGTTTGAGTGTTTATTGGTTTTCTTAAGAGAAATGGACTATTTTGCTCTCATGTGTAA
TGTTTTCTATTTAAATCTTTCTTAAATATACCAGCTGTTCTTTCCCTGAACTCTCCCCAGGTTCTA
GGACAAATTTAATAACATGTAATTCTCCTCAAATACTTTTGATGTCTCAGTGTTGGTGTTTTCTCCC
TAAACTAACATTAGGGCTGTGCCACGGGCATGACTTTATTTTGTGGGCTTTTTTTCCCTGCTTAA
GGAGAGGTGTCTTTTGGATATGAGCTATTTATTTGTGAAATGAAATTTGTCACCCAAATGATTCT
CTTATAAACTATTTGTAATGTCACTTATTCATTAGTGTTTGACATAATTTTAGAATATTTATTTGA
ATCAATCCTTTTATTACGAAAGACTTGAAGTTTGTGTCCATTCTTACAAGCCCTGGTCAGTCAAGTCC
CAATAAATGGTCAGCACAAAAAGATATTCTTGAAATTTGCTCTTTATTAAGGTATTACGTTGAGTTTG
CAACCAGATGGGAAAAATCACAAAAATGAGAAGGGGAGCAGATATCTTGTGAGGTCTGGATATTATCT
CCCTTTATAAACTTGGTGTAGGCCGTATATTTGAAAAATAAACATCTGGTTGAGGTTATTTCAATTTGGA
AACCTTCTTAGAGTCTGTGCTCAGTGTTGGTCCATCAGACGGTAGTACTGCCACAAGCTGAGAGCTT
GCTGGGCACGCAACCCAGAGTTTCTGAATCCTGTTCCGATTTTAACTTCCCAGGTGATTCACAGCCA
CATTAAGTTCAGGAGACCATAGGACCCCTGACCCATGAGAAGAAGTGATTGGTTTCGATTATAGACAT
CAAGATAAAGATTTTCTCAAAGGAAAGAGATTAATGTATTAACAACCTCTGAAAAATC
ACGCGTAAGCGGCCGCGCATCTAGATTGGAAGAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
  
```


[View online »](#)

Restriction Sites:	Sgfl-Mlul
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_006358.4</u>
Summary:	This gene encodes a peroxisomal membrane protein that belongs to the family of mitochondrial solute carriers. It is expressed in the liver, and is likely involved in transport. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Sep 2013]
Locus ID:	10478
MW:	46.4