

## Product datasheet for SC213249

### RFC3 (NM\_002915) Human 3' UTR Clone

#### Product data:

Product Type: 3' UTR Clones  
Product Name: RFC3 (NM\_002915) Human 3' UTR Clone  
Symbol: RFC3  
Synonyms: RFC38  
Mammalian Cell Selection: Neomycin  
Vector: pMirTarget (PS100062)  
ACCN: NM\_002915  
Insert Size: 1235 bp

Insert Sequence: >SC213249 3'UTR clone of NM\_002915  
The sequence shown below is from the reference sequence of NM\_002915. The complete sequence of this clone may contain minor differences, such as SNPs.  
**Blue**=Stop Codon **Red**=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
GAGGATGGATTGGAAGGCATGATGTTCTGACTTCTGTCTGAGTTATTCTTGCAAAGATTTCTCAGTATCAG
TATTTACATACAGCTTATATTTAAAAGAGCTGTGGGTAAATTAAGTGAAGTAAATCATGTCGTATTTGCG
TTTTTTTGGTAATAACTTCTCTGTGAAGTATAATCATCCTCTGAGTTAAATAATTGCTCCTATACTAT
TGAAGTATGTAGTTTTGTACATAACTTAGAGACTTTAGAGTCTAAGAAAATGATCTTAATTTACTTTAA
GCATTGGTTATTCAAGTATTCATTGTTGATCCTCCTATTCTCTCCGTCTAATCTCTCACCTGCTAAAG
GAGATTTACACATTAGAAAGCAAAGATTATTTTCATTTATCCAGATGACCATTTTCTGCCACAGGTAAC
ATGATTGTTTGACACACCATTATATTTAATTCTAGTTTCTCTCAATGAATAAATTGATTTTTGTAGGAA
ATGTAAGATTTCACTTCTGAAACATAATTATTGGTATGGACAAAATTGCAGATACCATTTCTGTTGAGGC
TGCAGATTTCCAACCTTTATTTCAAGTGGTTCAGATTAGTATTAGGTCGGTACTAAGAAAATAGCATGTT
TTCACATAATTAAGTACTTGAGACTCTTGAAGAAAATTCAGAATGAAGTCTGGAGAAAGGTATGTTAC
TGTAGTAATTAATCTTTGAACAGGTTTTGTGTTTTGTCTTAGCTCTGCCCTTTTTAATTAATATTTT
GGTTCATGGACCAAAGGTTTTACTTGACAAAATTTGTGTGACAGACTCCGAACAATTCCTTTACTACGAA
GTATAATTTATAAAAATAAATATACCCATTTTAAGGGTACAGTTTGATTTTTGACCAAGTAAACTATGA
TCCCAATCAAGGTATAGATGCCGTCACCCCAAAAAGTTCCCTCCATATCCCTTTGCAAGTCAAGTTCATCC
CTACCCTGGCCCAGATGATCACTGATCTTGTCTTATAGATGAGTTTTGCCAGTTCAAGAATTTAATGG
AATCAGATATTGTAAGCATTCTTGTGTAATACTTCATTCTCTCATTATTGAGATTCATCCATATTGT
TGAATGTTTCACTAGTTAATGTTTATTGTTCAATATTTTTGTATATACTTTAAAGCCTATTCACTTGC
TGATGGATCTTTGGTTTTGTTCCAGGTTTGGTTATTATGAATAAAGTTGCTGTGAATACTTA
ACGCGTAAGCGGCCGCGCATCTAGATTCAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
```



[View online »](#)

<b>Restriction Sites:</b>	Sgfl-Mlul
<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
<b>Components:</b>	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
<b>RefSeq:</b>	<a href="#">NM_002915.4</a>
<b>Summary:</b>	The elongation of primed DNA templates by DNA polymerase delta and DNA polymerase epsilon requires the accessory proteins proliferating cell nuclear antigen (PCNA) and replication factor C (RFC). RFC, also named activator 1, is a protein complex consisting of five distinct subunits of 140, 40, 38, 37, and 36 kDa. This gene encodes the 38 kDa subunit. This subunit is essential for the interaction between the 140 kDa subunit and the core complex that consists of the 36, 37, and 40 kDa subunits. Alternatively spliced transcript variants encoding distinct isoforms have been described. [provided by RefSeq, Jul 2008]
<b>Locus ID:</b>	5983
<b>MW:</b>	48.7