

Product datasheet for SC213256

MMRN1 (NM_007351) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	MMRN1 (NM_007351) Human 3' UTR Clone
Symbol:	MMRN1
Synonyms:	ECM; EMILIN4; GPIa*; MMRN
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_007351
Insert Size:	1226 bp
Insert Sequence:	>SC213256 3'UTR clone of NM_007351

The sequence shown below is from the reference sequence of NM_007351. The complete sequence of this clone may contain minor differences, such as SNPs.

Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAACGATCGCC
TTTAGTGGCTATTTATTATATCGTACAATAGTTAGTATGAAAAACAGACTATCACCTTTATTGAGAAAC
AGCCAGTGTTCATTTATCTTTGCTTGACATCTGCTCTGTTTTGGTTTTCTACAGGAAATGAAAT
CAACTGTTCATTTAATATGAGTAACTTGTATGTCTATTTTATAAAATTATTTGAATATTGTTAATG
TCTGAATATGAAAGAGTTCTTGATCCTAAAGAAATTTAGTGGCACAGAAAACAAAGTGAATTTGTTAGC
ATAATTATTCCTATTCTTATTTCTTCATTTAAGTCATTGCAATGGAAAGTAATATTATAAAATGGTAA
TTACAACATATTATCAGTCACAGTTTCTTTCCAATTAACACTTAACCTTTGTTATTCCCTGTATATA
AATATATAACACACATTTTCTAGATTACAAAATTTAAATAAAATTAACCAAAAAATGAAAATTGATTTTG
TAACTTTTATTTTACTCTTTACGTTGAGTTGATCAATTTCCATACTAAGATTTTCATTCAGAATCA
AAATTAAGAAAGTTGGACTGAAAATATGAAAAATGCTTAATTTCTCTCTCTATAATTCTCTAAT
TATAACATAGTAATTTACATGTAGTTGGACATGTACACTCAAGTCTAAGAATATATGAGTGGATCATTT
ACCGCCCCCGCCCCACAACATCTATAAGGGGCAAAAAGTCTTTTCTAATAAGTATTCTCTATGGTA
GTACCTACAGATCTGCCCTTCTTCTCTAAAGGTAAGTCATAATCTGTGTAATACTACAATTTATGGG
ATGCTCACTATGCCCTGTTTCTCTCTAAACAATTTACATGTAATGTCTCATTCCCTCACAATAACCCCTT
GTAAAGTGGGCATGATTACCATGATTTTATAGTTGAAGAACCTAAGACACAGAGACCAAGGCCCATGA
GCTCATAGGGCTGAGGCAGGATTTGGAATCAGGCCATGTCTTCTCCAGAGCCACATCCATCCTTTCTC
TATATTGCCTCCACAGATGTGCTAAAATTTATTTAACTAATCCTTTATCCTCTATTTGTGTGTCTCC
CATTTTTTATTATTACAATATTACTGTGGTGAACATGCTTAAAAATACATTCCTTGGATATCTGACAAC
GTGTTTCTGAAAAACAGATTTTCATAAGTAATAATAAAAAATAATAAAAAA
ACGCGTAAGCGGCCGCGCATCTAGATTCTGAAGAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
  
```


[View online »](#)

Restriction Sites:	Sgfl-Mlul
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_007351.3</u>
Summary:	Multimerin is a massive, soluble protein found in platelets and in the endothelium of blood vessels. It is comprised of subunits linked by interchain disulfide bonds to form large, variably sized homomultimers. Multimerin is a factor V/Va-binding protein and may function as a carrier protein for platelet factor V. It may also have functions as an extracellular matrix or adhesive protein. Recently, patients with an unusual autosomal-dominant bleeding disorder (factor V Quebec) were found to have a deficiency of platelet multimerin. [provided by RefSeq, Jul 2008]
Locus ID:	22915
MW:	47.8