

Product datasheet for **SC213273**

CHCHD7 (NM_024300) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	CHCHD7 (NM_024300) Human 3' UTR Clone
Vector:	pMirTarget (PS100062)
Symbol:	CHCHD7
Synonyms:	COX23
ACCN:	NM_024300
Insert Size:	1404 bp



[View online »](#)

Insert Sequence: >SC213273 3'UTR clone of NM_024300
The sequence shown below is from the reference sequence of NM_024300. The complete sequence of this clone may contain minor differences, such as SNPs.
Blue=Stop Codon Red=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
TTGAGAGCAGTGGGAAATATGCCCTATGAATGTTTGCATTAAGAGTGTATATAACTTAGAAGCAGA
TGAATATTTCTAATAAATGATTGCTGTAATATTTAAGACTGTACACCCCTCACCCAGACAGACCTTAAG
TTCTTCAAGTGGAGACAGTGAAGTCACCCCGTGTCTTTTTGCTTGTCTCAGTGCCATGCCGATGGTA
TGTGCTGTTGGCTGTGTGTGGCATGAGTTTGCATGACTTTCTGGAGGCATGGAGTTAGGTAAGGCTA
CATGAGAAATTGAGCTTTTCCACTGGGTTTTGAAAGAAGAGTATGATGTGAACAAGTAAAGACTGAATG
GGGCTGAGATGAAGGCAATGTTTCCAAGGAAAGGAAATGTTATGAGCAAGAGTGTGAGGCAAGAGAAGC
TGGAAACCACATTCAGAGAGTATCCTGTAGATTGCTCCACCTAGAATCTCAGGTGGGTGGAGCAGTGGTG
GGAGAAGACTGGAAAGGTAAGTTGAAGGTAAGGAATGTGTGGTGGCCTCAGATCCCAGGCTCATTCTT
CAAATCACTTCTTACTTCCCTCACTTATCTTTGTTTAAATAAGGTTAGTACATTCCTAGGGGCAATG
GGTTTTTCTAATAAATGACATAAAAAAGAACTGGTTGTCTTTCTTTAAAGAAAGTCTTAGTTACAGGA
TTCTGTGATTAATAATCCTTAAATATTGGGTTGCCTTCTGCAGCTGTAGTATCAGTTTTTGAACCACAGT
AATGGGAAGAAAGACAAATGGATTCCTTAGAAGTAAGATTTCTATTTGCAGGATGAGTTGGGCAGGGA
AAAGGGTCAGGGTTCATCAGGTGAACCTCAACACTGGGATGAGACTAGAATTCACCTTATGATATAAAC
ACAATACGATTATTCAATGTGGTACTGGGTAGACTGTGAAGCAGCCCTCCTTAGCCACTGCCTTT
TAGAATCGTTTGTATTCATGGTAGTTTATGAAGACATACTATTATTGAATGAAATTAATGTGTA
CTTGAAAACATTGCTTTTGTCCCTTCTCTTCATCTGGTCTTGGGTCAAGAACATGTTTTAATGGCTGC
CGACAATGAAGTCTGTCTGAGTCTAAAACCAAGCTCAGGTTTCTAAGCCACATGACCTTGATTGTTA
AATGTAGTTATATTTGAATAAAAAATGAATCATTTCAAAATCTTTGATTTGTATTGACTGCTATAA
AAATGATGGATTAATAAGGACCAGAAAAGAAATGAATGCCTACCTGTGTAAGAAGTATTTTTGTATTTT
TAAGCACACTTGTAAAGGGCTTGAATTTTAATGTCACCTCTTGAAGGTACATACCAATATCAACAT
TCATCTTGTAAGTGTCTCTGA
ACGCGTAAGCGGCCGCGCATCTAGATTCAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTTGATTCACCGCCGCTTCTATGAAAGG
```

Restriction Sites: SgfI-MluI

OTI Disclaimer: Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).

Components: The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.

RefSeq: [NM_024300.5](#)

Locus ID: 79145

MW: 53.9