

## Product datasheet for **SC213274**

### CHCHD7 (NM\_001011668) Human 3' UTR Clone

#### Product data:

Product Type:	3' UTR Clones
Product Name:	CHCHD7 (NM_001011668) Human 3' UTR Clone
Vector:	pMirTarget (PS100062)
Symbol:	CHCHD7
Synonyms:	COX23
ACCN:	NM_001011668
Insert Size:	1404 bp



[View online »](#)

<b>Insert Sequence:</b>	<p>&gt;SC213274 3'UTR clone of NM_001011668 The sequence shown below is from the reference sequence of NM_001011668. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site</p> <pre>GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC TTGAGAGCAGTGGGAAATATGCCCTATGAATGTTTGCATTAAGAGTGTATATAACTTAGAAGCAGA TGAATATTTCTAATAAATGATTGCTGTAATATTTAAGACTGTACACCCCTCACCCAGACAGACCTTAAG TTCTTCAAGTGGAGACAGTGAAGTCAACCCGTGCTTTTTGCTTGTCTCAGTGCCATGCCGATGGTA TGTGCTGTTGGCTGTGTGTGGCATGAGTTTGCATGACTTTCTGGAGGCATGGAGTTAGGTAAGGCTA CATGAGAAATTGAGCTTTTCCACTGGGTTTTGAAAGAAGAGTATGATGTGAACAAGTAAAGACTGAATG GGGCTGAGATGAAGGCAATGTTTCCAAGGAAAGGAAATGTTATGAGCAAGAGTGTGAGGCAAGAGAAGC TGGAAACCACATTCAGAGAGTATCCTGTAGATTGCTCCACCTAGAATCTCAGGTGGGTGGAGCAGTGGTG GGAGAAGACTGGAAAGGTAAGTTGAAGGTAAGGAATGTGTGGTGGCCTCAGATCCCAGGCTCATTCTT CAAATCACTTCTTACTTCCCTCATTATCTTTGTTTAAATAAGGTTAGTACATTCCTAGGGGCAATG GGTTTTTCTAATAAATGACATAAAAAAGAACTGGTTGTCTTTCTTTAAAGAAAGCTTAGTTACAGGA TTCTGTGATTAATAATCCTTAAATATTGGGTTGCCTTCTGCAGCTGTAGTATCAGTTTTTGAACCACAGT AATGGGAAGAAAGACAAATGGATTCCTTAGAAGTAAGATTTCTATTTGCAGGATGAGTTGGGCAGGGA AAAGGGTCAGGGTTCATCAGGTGAACCTCAACACTGGGATGAGACTAGAATTCACCTTATGATATAAAC ACAATACGATTATTCAATGTGGTACTGGGTTAGACTGTGAAGCAGCCCTCCTTAGCCACTGCCTTT TAGAATCGTTTGTATTTATCATGGTAGTTTTATGAAGACATACTATTATTGAATGAAATTTAATGTGA CTTGAAAACATTGCTTTTGTCCCTTCTCTTCATCTGGTCTTGGGTCAAGAACATGTTTTAATGGCTGC CGACAATGAAGTCTGTCTGAGTCTAAAACCAAGCTCAGGTTTCTAAGCCACATGACCTTGATTGTTA AATGTAGTTATATTTGAATAAAAAATGAATCATTTCAAAATCTTTGATTTGTATTGACTGCTATAA AAATGATGGATTAATAAGGACCAGAAAAGAAATGAATGCCTACCTGTGTAAGAAGTATTTTTGTATTTT TAAGCACACTTGTAAAGGGCTTGGAAATTTAATGTCACCTCTTGAAGGTACATACCAATATCAACAT TCATCTTGTAAGTGTCTCTGA ACGCGTAAGCGGCCGCGCATCTAGATTCAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA CGAGATTTGATTCACCGCCGCTTCTATGAAAGG</pre>
<b>Restriction Sites:</b>	SgfI-MluI
<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
<b>Components:</b>	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
<b>RefSeq:</b>	<a href="#">NM_001011668.3</a>
<b>Locus ID:</b>	79145
<b>MW:</b>	53.9