

Product datasheet for **SC213284**

ZMAT1 (NM_001011657) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	ZMAT1 (NM_001011657) Human 3' UTR Clone
Symbol:	ZMAT1
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_001011657
Insert Size:	1250 bp
Insert Sequence:	>SC213284 3'UTR clone of NM_001011657 The sequence shown below is from the reference sequence of NM_001011657. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAACGATCGCC
CTTTGGGATGAGTCTATTCTTGGATTTTAATGTTTAGTTTTGTTTACCCAAGGTTGAATTGAAAAAA
AAAACAGTCAATATGGATTTAGAAAAAGGAACACCTGATGAAGAAAGGAGAGGTAGATACAGTCAGTG
TCACTTCAGGACACTTAGGTTTTTTTTGTATAAAATTTAAATTGAATTTAAAGAAGGAAAAAAGC
CCAAACTTAACCTCTGAGAAAAGAACATAAGAAGTCAAGGAGAACATAAGAGAAAAGGAAACCTGTTAC
AGAAAAGACAAGAATCTGTGTTTTGGAATGAGTCTATTCTTGGGTATTGAACTTTTAGTTTTGTTTGCC
CAAGGATTAATTGAGGAATCAGCTAAGAAATGGACTTTAGACAAAAGCAAGAGGATCAGATGAAGAA
AAGGAGAGGTAGATACAGTCAGTGTCACTTCAGGAAAGCTATTTAAAAAACTTGAAATTTAATTGAAA
GAAGAAACAACAACAAAAAGCCTAAACCTAGCCTCTGAACAACATAACATGAGAACACAAGAAGCTTA
AGAGAAAAAGAAACCTACTCAAGAAAAGACAGAAGAGACAGTGATTTGGGATGAGTCTACTCTAGGATT
TTCACTTTTTAGTTTTGTTCTTCAAAGTTGAAGGAAAAAAGTTTGGTTTTATAAAATTCATGTTAT
TGTAATTTTTCTAGGTGGATGGCTACTTTAATCTCTAAAAAGCCAAGTGAAGTAAAGTATTTCAGTAT
GCCTTTTCTCAAGTTACTTTCTTCTATTTCTTAAAAAGAAAAAAATTATTAATGTTTCTCACAT
ATCTCACATATAATGTAATTTCCCTAAATGAAGTTGTCTCTACTTCTGCTCATCAAATTGCTGTGATAG
TGAATTATTTATTCATGGGAGATAATTTATTTTAAAGGACAGAATTACCAAGCGTTACAAAATCAGTTC
TTTCCTTGGTTTTGTGTTAGTGTGGTGGTATTTTATTGTTGTTTTCTGTGTTTATGTGTCTCAGCTT
TCTCCAAGGAATATGTATGAAATAACTTAACTGATTTTTCTTGTAAATCTAATTTGCAGTGTATT
TTTGCAATTTCTAGTTCTGAAAGTGAAATGAAACAGTCTATAATAAAGTATAGATATATAGTTTT
AAAACGGTCTCAAAAAGTACTGATATAAGGTCAGTCTATATTCTGGAAATGTTTATTTAAAGTGTTTT
AATTTCTA
ACGCGTAAGCGGCCGCGCATCTAGATTGAAGAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTTCGATTCCACCGCCGCTTCTATGAAAGG
  
```


[View online »](#)

Restriction Sites:	Sgfl-Mlul
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_001011657.4</u>
Summary:	This gene encodes a protein containing Cys2-His2 (C2H2)-type zinc fingers, which are similar to those found in the nuclear matrix protein matrin 3. Alternatively spliced transcript variants have been found for this gene. [provided by RefSeq, Jan 2012]
Locus ID:	84460
MW:	48.9