

## Product datasheet for SC213328

### THNSL1 (NM\_024838) Human 3' UTR Clone

#### Product data:

|                           |                                       |
|---------------------------|---------------------------------------|
| Product Type:             | 3' UTR Clones                         |
| Product Name:             | THNSL1 (NM_024838) Human 3' UTR Clone |
| Symbol:                   | THNSL1                                |
| Synonyms:                 | TSH1                                  |
| Mammalian Cell Selection: | Neomycin                              |
| Vector:                   | pMirTarget (PS100062)                 |
| ACCN:                     | NM_024838                             |
| Insert Size:              | 1239 bp                               |

**Insert Sequence:** >SC213328 3'UTR clone of NM\_024838  
 The sequence shown below is from the reference sequence of NM\_024838. The complete sequence of this clone may contain minor differences, such as SNPs.  
 Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAACGATCGC
GAACAACCTTGCCAAATCAATTCATATGAAAGCTTTCAGAGTAAATTTTTTTTCTAGCTATAAGCAT
GCAATAATAAATCTCAAACTGATTTGGAGTACAGTAGCATTTTGTCTTTATGTAAATATCTCTATA
TCTGTTTGAATTTCAAAGTCTGATCTCTAGGGCTGACATGAGAACTCCATGTCACCTCCTCAGATCT
TTAATCTGGAAGTGACAAAAGGTAACACAGTGCACGGACCTTTGAGCTATCATACTTTTGCCTTCTGCT
CATGAGGGGTGTATCAATTTCCACTCCCACACCTCACTGTTATGTGGACCAAAATGTCTGGTATACTA
TTTGGCGATTAAATATTTAAGCCAGTTTTCAGGTAACATCTGTAACAGTGAATACTCTGTTGAT
TAGAAAGTTAATTTACCACTCAAATGGCTAAGTTGAATGGAGGAAGTAGTAAACTCTTCTATTAGGAT
TTGGGCCAGTGTTTTTGTGTTTTTGGTGCACTTTCAGCTCACTGCAACCTCCGCTCCCAGGTTT
AAGTGATTCTCCTGCCTCAGCTCCCAGGTAGCTGGGATTACAGGCATGCAACACCACACATGGCTAAT
TTTTTATTTTAAATAGAGATGAGGTTTCGCCATGTTTTCAGGCTGATCTTGAACCTCTGACCTCAGG
TGATCCATCCACCTCGGCTCCCAAAGTGCTAGGATTACAGGCGTGAGCCTCTGCACCTGGCTAGGGCC
AGTGTTTTTAACCCTCCTCAAATTCATTTTGGACTACAATAATGTGTATCTCACACTGTCCAAATTA
AAGATGTACTCTGTCTTCACTCTCCTTTGGATGTGTACGGTTAGGGTGACGATGTTTGAATCAACA
CATAACTGAAATCCAGGATGAAACAGTAAGCAGAACTACTTAGAAGTGAATACAGAAATGATCTATTTT
AAGAGGAAGGTTTTCCAAATTAAGAGTACTTGAGTAGTCTCAATAGGAGTGATTTGTAGACAGCAGT
TTCCCTATATTATTTGGAGTCAATTTACGTAATTTGTTGTTTATTTCTTCAGTGAACATTGTCTTTT
TCTGGCAGCTGAACTGCACAACTGATACACATATTTAATTTGTATTCTTTGTAGTAATACATTTT
AACCATGATTTATGTTTGTGAATTTACTGGATGTTTTCTTCTGAATTAAGATTGTCAAGACCAA
ACGCGTAAGCGGCCGCGCATCTAGATTCTGAAGAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTTGATTCCACCGCCGCTTCTATGAAAGG
  
```


[View online »](#)

|                           |                                                                                                                                                                                                                                                                                            |
|---------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Restriction Sites:</b> | Sgfl-Mlul                                                                                                                                                                                                                                                                                  |
| <b>OTI Disclaimer:</b>    | Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs). |
| <b>Components:</b>        | The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.                                                                                             |
| <b>RefSeq:</b>            | <u><a href="#">NM_024838.5</a></u>                                                                                                                                                                                                                                                         |
| <b>Locus ID:</b>          | 79896                                                                                                                                                                                                                                                                                      |
| <b>MW:</b>                | 47.6                                                                                                                                                                                                                                                                                       |