

Product datasheet for SC213335

RBM41 (NM_018301) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	RBM41 (NM_018301) Human 3' UTR Clone
Symbol:	RBM41
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_018301
Insert Size:	1236 bp
Insert Sequence:	>SC213335 3' UTR clone of NM_018301

The sequence shown below is from the reference sequence of NM_018301. The complete sequence of this clone may contain minor differences, such as SNPs. **Red**=Cloning site
Blue=Stop Codon

CAATTGGCAGAGCTCAGAATTCAA**GCGATCGC**

GCTACAGGTAGCACTACAGAAATTAGTGGTAGCT**AG**AATATATATAAATTGTGGGTCCTGGATGGTCTTT
 CTTGAATTGGAATCTTGGATCTAGGATTTATGTATAGACACCAAGTTTGGTGGAGTTGAAGAGGAAAAAC
 TAACAGAGGGAGGGGAAAAAGTAATTTTCAGCTTTTGGAGAGAACTGGTAGAAAACATTTCTATTATAT
 AATTGAAGACTATATAGGACAGATGATGCTGAGTATGAAGGGTAAAGTGCAAAGAGTATTATTTTCTC
 AGGATGAATTACAACGGATTATTCCTAAGGATTTTAAAAATCTTTATAATTTTCTGGGGTCCTCAATA
 TTTGTCTCTAGAGATAACATTCTTTCTTAGACTGCTTCGTATACACAGAAATAAAAACTAATTAATGT
 TGAAATATTTCTTACTAGAGCCCATTTGCCATATTTTCATCCACCTAGTAAATAAAGAGTAGGACCCCC
 AGGCCTAAGCATCAACAAATGCTATTAAGTCAGAACTATCAGCAGTAATAGTAGGTAGTTTCATAAGTCA
 GTTGACTTATTTAATGAATGCTTCCTCCACAAGTTATCAGTCTTCAAAAAGCTATATATTTTCTAGTGA
 AATATTTTAGATCAAAATGTTGGTGTCTCACAGCCATTTAAGCTAGCAACCATTGAGTATCATTGAATT
 CATAGCAATAACATTGATTAAGTGTCTAATATATGTTAAATATTAGGCATATGAATTTTGAAGAGGAA
 CATACAAGCAAAGAAATGCCAAACGGAAATGGAAATAAAGCGTTCTCTAGAAAAGAGGATGCCAGTCACAA
 ATCCCCAGTTGTTTGTATGTTGTCTGGAAAATAAAAAACAATATACAGAGAGGAAGTAGTACAAAC
 TGGAGTCAGACCTGTACAGCCACACTTAACATTCTATGACCTTGGCAAGTCACTTAATCTGAGCCTCACA
 TTTTTCATCTATAAAGTGAGGATAAATGCTACTTACCTCACTGGGTTGTGAAGATGAAGTGAAGATATT
 ATCCAGGCTGAAGAACAGAAAGGAAAACAATGAATAAAATTGAACAAGAGTATAGAACCTGTGTGACACC
 ATTAAGCATTCTACATATTTCTAATGGGAGTACAGGAAAAAGGAAAAGAAATGAAGGAGCAGAAAAAT
 ATTCTAAGAAATAATGGCTATAAATACCCCCGATTTGATGAGAGAC

ACGCGTAAGCGGCCGCGCATCTAGATTCAAGAAAATGACCG



[View online »](#)

Restriction Sites:	Sgfi-MluI
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	NM_018301.3
Summary:	May bind RNA.[UniProtKB/Swiss-Prot Function]
Locus ID:	55285