

## Product datasheet for SC213343

### PHR1 (PLEKHB1) (NM\_001130036) Human 3' UTR Clone

#### Product data:

**Product Type:** 3' UTR Clones  
**Product Name:** PHR1 (PLEKHB1) (NM\_001130036) Human 3' UTR Clone  
**Vector:** pMirTarget (PS100062)  
**Symbol:** PLEKHB1  
**Synonyms:** KPL1; PHR1; PHRET1  
**ACCN:** NM\_001130036  
**Insert Size:** 1247 bp  
**Insert Sequence:** >SC213343 3'UTR clone of NM\_001130036  
 The sequence shown below is from the reference sequence of NM\_001130036. The complete sequence of this clone may contain minor differences, such as SNPs.  
 Blue=Stop Codon Red=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
TCGCTCATGTGGTCGCCCTGCTGGTTCTGAGGCCCTGGGACTCGGAGCACTGACCCTGCGCTTGATTG
CTAGACTCCTCTTCCTCCTGGACCCATCCTCTACCATCCAAGCCCTGTCCCCTTTGGCCCTATCCTC
TCCATTAGCTCCTTCCGGGTTTGGACCATTCCCCCACTCCCTACCCTTAATCCCCACATGGGAAGAAG
CTATCATCACAGGTACAACATCGCTTGAAGTCTTACATCTACCACTAGACACCCCAAAATCTGTTA
TAGACATTTATGGATACATTTCTCTAAACACAACAGGGCACAGCAAATACGACTTCATTTGGCTTCGA
GTTCCCCAGGCGCTGTAGACACAACATGAATCGGGCTCTCTGCTCTCCTTAGGGAGCTCGAGTCTG
GTGGGGAGAACAGGAGTAAACAAGGACTTGACAAAGCTGAAGAGTTATCAGTCCTTTGACAAGGACAGG
TGGGGCAGGGAGCAAGACAGGTAGGCTGGAAGAACAGTTATTGGCAAGTATGCAGAGCCGTGAACGTCA
TGGCATGTCCAAGGAATTAATGGGAGTTCAATTTGGGCTGGGGTGGAGGCTGGGATCAGACCGTGGTGG
GCCTTCAAGCTAAGGAGCTTCTAGGTGAAAGGGGAGATGTGAGCCTTCTCTGGAGGGAAGTTTCATGA
TTGCATCTATAATGAATATATTGCTGTTTTGTGAATACTGACACATGTCCATACCTAAAACACTCCTG
AGTCAAGTCCCATCCTTCCCACAACAGCTTCTGGCTGGTACCCATGATAACAATTGAGCTGAACCTG
GGGACCCCTGGTTGGGGAACAGGTGAGTTCTATTTGAGACTTCCAGCCCTAGAAAGCTGCCTCCGTTCCA
GAAATGCCTCTCACACCAGGAGCTCGGCCCTCTTTTGTAGCTGTGACTGTACCCTCTCAGGCTTTGT
CTCATCCTTATTCTGAATAAGATGGCAGTGTTCTCCTCTGGGCGCTGATCCACCTCTACACCAGCCCA
GGAAGCCCCATCTGTGCCTGCCCTCAGGTGGTCCACCAGTCTCCCCCTTTGGTTCCCTTCCAGTCTCTT
CCCCCTTCTATCCCAATCACCAATAGAAATGCTAACATCCCTGCCTGGTAGCCAGACTAGCCCACTAA
AGCTCCCCTGTAATGGGGGCTCCATTAGTTCTGCTGCCGAGACTAATAAAGATTTGGTTGGCTCTAGC
AGTAA
ACGCGTAAAGCGGCCGCGCATCTAGATTGAAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCCTTCTATGAAAGG
```

**Restriction Sites:** SgfI-MluI



[View online >](#)

<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
<b>Components:</b>	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
<b>RefSeq:</b>	<u><a href="#">NM_001130036.2</a></u>
<b>Locus ID:</b>	58473
<b>MW:</b>	46.1