

Product datasheet for **SC213344**

PHR1 (PLEKHB1) (NM_001130035) Human 3' UTR Clone

Product data:

Product Type: 3' UTR Clones
Product Name: PHR1 (PLEKHB1) (NM_001130035) Human 3' UTR Clone
Vector: pMirTarget (PS100062)
Symbol: PLEKHB1
Synonyms: KPL1; PHR1; PHRET1
ACCN: NM_001130035
Insert Size: 1247 bp
Insert Sequence: >SC213344 3'UTR clone of NM_001130035
 The sequence shown below is from the reference sequence of NM_001130035. The complete sequence of this clone may contain minor differences, such as SNPs.
 Blue=Stop Codon Red=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
TCGCTCATGTGGTCGCCCTGCTGGTTCTGAGGCCCTGGGACTCGGAGCACTGACCCTGCGCTTGATTG
CTAGACTCCTCTTCCTCCTGGACCCATCCTCTACCATCCAAGCCCTGTCCCACTTTGGCCCTATCCTC
TCCATTAGCTCCTTCCGGGTTTGGACCATTCCCCCACTCCCTACCCTTAATCCCCACATGGGAAGAAG
CTATCATCACAGGTACAACATCGCTTGAAGTCTTACATCTACCACTAGACACCCCAAAATCTGTTA
TAGACATTTATGGATACATTTCTCTAAACACAACAGGGCACAGCAAATACGACTTCATTTGGCTTCGA
GTTCCCCAGGCGCTGTAGACACAACATGAATCGGGCTCTCTGCTCTCCTTAGGGAGCTCGAGTCTG
GTGGGGAGAACAGGAGTAAACAAGGACTTGACAAAGCTGAAGAGTTATCAGTCCTTTGACAAGGACAGG
TGGGGCAGGGAGCAAGACAGGTAGGCTGGAAGAACAGTTATTGGCAAGTATGCAGAGCCGTGAACGTCA
TGGCATGTCCAAGGAATTAATGGGAGTTCAATTTGGGCTGGGGTGGAGGCTGGGATCAGACCGTGGTGG
GCCTTCAAGCTAAGGAGCTTCTAGGTGAAAGGGGAGATGTGAGCCTTCTCTGGAGGGAAGTTTCATGA
TTGCATCTATAATGAATATATTGCTGTTTTGTGAATACTGACACATGTCCATACCTAAAACACTCCTG
AGTCAAGTCCCATCCTTCCCACAACAGCTTCTGGCTGGTACCCATGATAACAATTGAGCTGAACCTG
GGGACCCCTGGTTGGGAACAGGTGAGTTCTATTTGAGACTTCCAGCCCTAGAAAGCTGCCTCCGTTCCA
GAAATGCCTCTCACACCAGGAGCTCGGCCCTCTTTTGTAGCTGTGACTGTACCCTCTCAGGCTTTGT
CTCATCCTTATTCTGAATAAGATGGCAGTGTTCTCCTCTGGGCCTGATCCACCTCTACACCAGCCCA
GGAAGCCCCATCTGTGCCTGCCCTCAGGTGGTCCACCAGTCTCCCCCTTTGGTCCCTTCCAGTCTCTT
CCCCCTTCTATCCCAATCACCAATAGAAATGCTAACATCCCTGCCTGGTAGCCAGACTAGCCCACTAA
AGCTCCCCTGTAATGGGGGCTCCATTAGTTCTGCTGCCGAGACTAATAAAGATTTGGTTGGCTCTAGC
AGTAA
ACGCGTAAAGCGGCCGCGCATCTAGATTGAAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCCTTCTATGAAAGG
```

Restriction Sites: SgfI-MluI



[View online >](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_001130035.2</u>
Locus ID:	58473
MW:	46.1