

Product datasheet for **SC213371**

TMEM140 (NM_018295) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	TMEM140 (NM_018295) Human 3' UTR Clone
Symbol:	TMEM140
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_018295
Insert Size:	1246 bp
Insert Sequence:	>SC213371 3'UTR clone of NM_018295 The sequence shown below is from the reference sequence of NM_018295. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAACGATCGCC
AGGGCTGAGAGCAAGCTTGAGAGCTGCTAAAGGCTTACGTGATTGCAAGGGTTCAAGTCCAACCATGGT
CAGAGGTGGCAGATCTGCTCAGCCATCTCATTTACAGCTAACGCTGATCTCCAGCTCCAGCGATGGAA
CCCACTACAGAGGAGGTGGGGCCCTGTGTCAAAGAGGCCGAGGGGCAGCAAGGGCAGCCAGGGCACCT
GTGACTTCTTAGTACAAGATTGTCTGCTTCAAGGCTCCAAGGCTCCCAAAGACTCCCTAAACCATG
CAGCTCATTGTACACCAATTCTGCTTTAATTAAATGGATCTGAGCAAATCTTCTCTAGCTTCAGGAG
GGTGGGGAGGGAGTGATTGCTGTCATGGGGCCAGACTTCCAGGCTGATTGCCAAATGCCAAATGAAA
CCTAGCAAAGAACTTACGGCAACAAACGAGGACATTAAGAGAGCGAGCACCTCAGTGTCTCTGGGGACA
TGGTTAAGGAGCTTCCACTCAGCCACCATAGTGAGTGGGCCGCCATAAGCCATCACTGGAATCCAAC
CCCAGAGGTCCAGGAGTGATCTCTGAGTGACTCAACAAAGACAGGACACATGGGGTACAAAGACAAGGC
TTGACTGCTTCAAAGCTTCCCTGGACCTGAAGCCAGACAGGGCAGAGGCGTCCGCTGACAAATCACTCC
CATGATGAGACCCTGGAGGACTCCAAATCCTCGCTGTGAACAGGACTGGACGGCTGCACACAAACAAAC
GCTGCCACCCTCACTTCCCAACCCAGAACTTGAAAAGACATTAGCACAATTACGATTGGGGAATTG
TGTGTATTTTCTAGCACTTGTGTATTGGAAAACCTGTATGGCAGTGATTTATTATATATTCCTGTCCA
AAGCCACACTGAAAACAGAGGCAGAGACATGTAATCTGCTGTGATCTCTTGTCTCAGTGTCTCTCTG
GGCTCCTGTCCCTCTTGCTTTATAGCTAGCTGCCCGGGGACCAAGGTACAGGTGAAAGCAAGGTAGCAG
CTTGCGGGAGGAGGCCTGTCTGGCTTACCAGTCTATACACTGTGGCCTCAACCTCCCAGACAGGGCAGA
GAACTGTGGGCAGCTCGTTTGCTTTCTAGGCTGGCTGGAGAGGTGGGAGCTCATTGATAGACTCATGAT
GGAAACTATTTTGAACAGGCTTCTCCTTCAGGAGAGATCATGCGGACTAACTGTAGCAATTCAG
TGCA
ACGCGTAAGCGGCCGCGCATCTAGATTGAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
  
```


[View online »](#)

Restriction Sites:	Sgfl-Mlul
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_018295.5</u>
Locus ID:	55281
MW:	45.7