

Product datasheet for **SC213372**

PRKAR1B (NM_001164759) Human 3' UTR Clone

Product data:

Product Type: 3' UTR Clones
Product Name: PRKAR1B (NM_001164759) Human 3' UTR Clone
Vector: pMirTarget (PS100062)
Symbol: PRKAR1B
Synonyms: PRKAR1
ACCN: NM_001164759
Insert Size: 1263 bp
Insert Sequence: >SC213372 3'UTR clone of NM_001164759
 The sequence shown below is from the reference sequence of NM_001164759. The complete sequence of this clone may contain minor differences, such as SNPs.
 Blue=Stop Codon Red=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAACGCATCGCC
TACAACAGCTTCATCTCCCTCACCGTCTAGGCACAGCTCCCGCCCTGCAGCCCCAGTCCCCAGTGTGG
TGGCCGTGCCCTGCTCGTCTGTGTGCGGGGCCCCGGGAGCCGCTGTGTGAGGTGTGGGCCGGTGGGGCT
GGGTCCCAGGAGCGTGTGAGACTGCCCTTCCCCGACTCACTTTTGAATAAATGATCACCTTGTGCA
TTTCCAAATCAAAGGACAAGCGGACAAAATGCATCCCAAGATCAAGGAAGGGACAGGCCAGCTTCTCTCC
CCACACGCTCCCCGGCTGCCTCTGTGGCTTCTCCTGGGGGGCCACCCACCCCTGCCAGTCTCTGTG
GAGATGCTTGTGAGGATCGGTCTCCCCAGAACAGGCCAGGACGTTGCCCTGGGGCTGTGACCCTGT
GAGGTCAGGTCAGGTCAGATTGAGGTCTGAGTGTGGCAAGTGTGTCAAAGGGGCTGCCCCAGGGGA
GATGAGGCTGAGAGCAGGAGTTGAGGCCAAGAAGTCAAGGCCCTCCCGCAAATGTGTACCCCTGCC
CGCGCCACTGCACCCCGCCGACCCCACTCCCCGGGGCCCTGCTGCGGATGCCGGAGTGGGAGAG
TCTCTGAGCTGTGAGATTGATCTTGCCCTAATTGGAGAGGAAGCCGGGGCCCAAGACACAGGGGCTC
CTGCCCTGGGAGCCAGGGCCGCGCCGAGGTAGACCCAGTAGGGGGGGCCGGGCTCGAAGTTCCTTT
GGGAGGGGCTGGCGGACTCCAGCAGGCCGTCTCACCTTTCTTAGAAAGTCCACCCAGGGCAAGTTGA
TGTTGGGGGAAAGCAGAAGTCAAGCCAGCCGCGGGCCCCACAGCCCCGACCCGCTCAACCGACTTGTC
CCTTAAATGTGTCTTGGATCCCGCAGTGATGACGTGAGCCAGCCAGGCCGAAAGGGTGAGGCCAGTGC
AGAGAAGCTTCCAGGGGATTCTGGGTCCCCGAAAGACAAGCAGGTCATTGCGATTCACCTGATG
TTGCTCCTGTCCCGTGCCTCCCGGGGAGGCTGTCTTGGTCCGAGGCTGCTCGTGCAGCCCTCCCT
GCTGGCGGTTACGATGCGTGGGGTCCCCCTCCCCACCCAGCCCCGGCCACGTGCGCGTGTCCCGCTGT
GCCTGTGACGTCCTGTGACCTGTGAGCCATCCCCCCTTATCTCTGCTCTGAATACTGCATGAAGC
ATTAACGTGCAATGAAGTAC
ACGCGTAAGCGGCCGCGCATCTAGATTCAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
```

Restriction Sites: Sgfl-MluI



[View online >](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	NM_001164759.1
Summary:	The protein encoded by this gene is a regulatory subunit of cyclic AMP-dependent protein kinase A (PKA), which is involved in the signaling pathway of the second messenger cAMP. Two regulatory and two catalytic subunits form the PKA holoenzyme, disbands after cAMP binding. The holoenzyme is involved in many cellular events, including ion transport, metabolism, and transcription. Several transcript variants encoding the same protein have been found for this gene. [provided by RefSeq, Aug 2015]
Locus ID:	5575
MW:	44.6