

Product datasheet for **SC213385**

WDR47 (NM_014969) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	WDR47 (NM_014969) Human 3' UTR Clone
Symbol:	WDR47
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_014969
Insert Size:	1246 bp
Insert Sequence:	>SC213385 3'UTR clone of NM_014969

The sequence shown below is from the reference sequence of NM_014969. The complete sequence of this clone may contain minor differences, such as SNPs.

Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAACGATCGCC
ACTGTCACCCCTCTGGACTTACAATGGGAGAGCACACCGCATGTCAGTCTATGCAGCAAAAGCACAGAG
ACTTAAGACTACTGAGTTGTGAAAATTACAAATCTGAAGAACATAGTGTCAGGAAAGTGGTTTAGCAC
GAAGAGGCCCTTATTACCATGTATCCCACTGATAGGAGGTGTTGGTGGTGTATTCCGCAGTGCTTT
CAGTCTTCCATGTGAGCTCGTGCTGCTGTGACCTGCTATATGTAGTCTCGTTGCCAAAGTCTGCAGAAG
AGCTCTTCAGTTGTTGGTGTGCACTCCAAGTCAGGATGGACAATGTGTTTACGGTTTAGTATTTCAATGC
ATTCTTGGTCTTTGCCTAAATAACAGTTTTATATGCACATTGAAATGGAATTATACTTCAACTATATT
ATTAAATGTAATGCAACCAAGTTCCTCCAGATTAAACTTCCAGGTGTTTCAAGATTACTTTTGCTCTT
CTCAGATCCCATATTGTATTATCACTTGTCTTAGAGGTGAGAATCCATAATATATGCACTCAAA
AGTTACATGGTTGCTTTCACTTAAGGATCATTGTGGAGTTTAAAGATGAATGAAAACTGCTTCTAGT
TTACTACATGGTATAGGCCCTTTTTCTTAAACCCAGGGATATGATTATTTTGCATATAATTTGTTT
CAGGCTAAAAGGTAAATGTGTTTGCTTCAGAACTTGTTAACTTCAGTTTTTTGAATGCAACAGGATAC
CTCCCTTCCAACTGAAGTGTAGAAGCAGAGCAGCAGGTTATGTGATGCAACACTTGATGGTACAGT
AAATTTACTGGCATTCTCTCTTAAAAATTAATCCTTGACATAGACCATAGCATGGCTTGAAATGC
TATGTCTGCATGATAATTTAAATGGAAGATTTAACTTTGCACTCCAAAGCTTATTTGGATTTTTTT
CTTGCACTGTTTTGTGTAATGCAGAATAATGATTTTATTTCTACAGCTTTGTAGATTCTAACATTTATG
TATCTTTATTTTCATATTGTACAGTAATTTTACTTTAAATTATTTAAATAGGCTATTTTATTTATTTCA
AATGCAGTTGTATTAGTTCTCATTATTGAAGTGTCTGTGCACTGTATGTAGCAAGCATTTTTCATCTGT
TGATATACAAGTGGAAGGGTATTAGAAGTGAAGTGTCTATTATTTCAATAAAGACCTCTTGACATTT
AAAA
ACGCGTAAGCGGCCGCGCATCTAGATTGAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG

```



[View online »](#)

Restriction Sites:	Sgfl-Mlul
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_014969.6</u>
Locus ID:	22911
MW:	48