

Product datasheet for **SC213387**

WDR47 (NM_001142550) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	WDR47 (NM_001142550) Human 3' UTR Clone
Symbol:	WDR47
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_001142550
Insert Size:	1246 bp
Insert Sequence:	>SC213387 3'UTR clone of NM_001142550 The sequence shown below is from the reference sequence of NM_001142550. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAACGATCGCC
ACTGTCACCCCTCTGGACTTACAATGGGTAGAGCACACCCGATGTCAGTCTATGCAGCAAAAGCACAGAG
ACTTAAGACTACTGAGTTGTGAAAATTACAAATCTGAAGAACATAGTGTCAGGAAAGTGGTTTAGCAC
GAAGAGGCCCTTATTACCATGTATCCCACTGATAGGAGGTGTTGGTGGTGTATTCCGCAGTGCTTT
CAGTCTTCCATGTGAGCTCGTGCTGCTGTGACCTGCTATATGTAGTCTCGTTGCCAAAGTCTGCAGAAG
AGCTCTTCAGTTGTTGGTGTGCACTCCAAGTCAGGATGGACAATGTGTTTACGGTTTAGTATTTCAATGC
ATTCCTTGGTCTTTGCCTAAATAACAGTTTTATATGCACATTGAAATGGAATTATACTTCAACTATATT
ATTAAATGTAATGCAACCAAGTTCCTCCAGATTAAACTTCCAGGTGTTTCAAGATTACTTTTGCTCTT
CTCAGATCCCATATTGTATTATCACTTGTCTTCTAGAGGTGAGAATCCATAATATATGCACTCAAA
AGTTACATGGTTGCTTTCACTTAAGGATCATTGTGGAGTTTAAAGATGAATGAAAACTGCTTCTAGT
TTACTACATGGTATAGGCCCTTTTTCTTAAACCCAGGGATATGATTATTTTGTATATAATTTGTTT
CAGGCTAAAAGGTAAATGTGTTTGCTTCAGAACTTGTTAACTTCAGTTTTTTGAATGCAACAGGATAC
CTCCCTTCCAACTGAAGTGTAGAAGCAGAGCAGCAGGTTATGTGATGCAACACTTGATGGTACAGT
AAATTTACTGGCATTTTTCTCCTTAAAAATTAATCCTTGACATAGACCATAGCATGGCTTGAAATGC
TATGTCTGCATGATAATTTAAATGGAAGATTTAACTTTGCACTCCAAAAGCTTATTTGGATTTTTTT
CTTGCACTGTTTTGTGTAATGCAGAATAATGATTTTATTTCTACAGCTTTGTAGATTCTAACATTTATG
TATCTTTATTTTCATATTGTACAGTAATTTTACTTTAAATTATTTAAATAGGCTATTTTATTTATTTCA
AATGCAGTTGTATTAGTTCTCATTATTGAAGTGTCTGTGCACTGTATGTAGCAAGCATTTTTCATCTGT
TGTATACAAGTGAAAGGGTATTAGAAGTGAAGTGTCTATTATTTCAATAAAGACCTCTTGACATTT
AAAA
ACGCGTAAGCGGCCGCGCATCTAGATTCGAAGAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTTTCGATTCCACCGCCGCTTCTATGAAAGG

```



[View online »](#)

Restriction Sites:	Sgfl-Mlul
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_001142550.2</u>
Locus ID:	22911
MW:	48