

Product datasheet for **SC213394**

CD45 (PTPRC) (NM_080921) Human 3' UTR Clone

Product data:

Product Type: 3' UTR Clones
Product Name: CD45 (PTPRC) (NM_080921) Human 3' UTR Clone
Vector: pMirTarget (PS100062)
Symbol: PTPRC
Synonyms: B220; CD45; CD45R; GP180; L-CA; LCA; LY5; T200
ACCN: NM_080921
Insert Size: 1325 bp
Insert Sequence: >SC213394 3'UTR clone of NM_080921
 The sequence shown below is from the reference sequence of NM_080921. The complete sequence of this clone may contain minor differences, such as SNPs.
 Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
GCAAGTCCAGCTTTAAATCAAGGTTCAAGGAAAAACATAAATGAGGAACTCCAACCTCCTGTTAG
CTGTTATTCTATTTTTGTAGAAGTAGGAAGTAAAAATAGGTATACAGTGGATTAATTAATGCAGCGA
ACCAATATTTGTAGAAGGTTATATTTTACTACTGTGAAAAAATATTTAAGATAGTTTTGCCAGAACAG
TTTGACAGAGCTATGCTTATTTTAAAAATTTATCTCTTATTCAGTAAAAACAACCTCTTTGTAATCG
TTATGTGTGTATATGTATGTGTGTATGGGTGTGTTTGTGTGAGAGACAGAGAAAGAGAGAGAATTCT
TTCAAGTGAATCTAAAAGCTTTTGCTTTTCCTTTGTTTTATGAAGAAAAATACATTTTATATTAGAA
GTGTTAACTTAGCTTGAAGGATCTGTTTTTAAAAATCATAAACTGTGTCAGACTCAATAAAATCATGT
ACATTTCTGAAATGACCTCAAGATGTCCTCCTTGTCTACTCATATATCTATCTTATATAGTTTACT
ATTTTACTTCTAGAGATAGTACATAAAGGTGGTATGTGTGTGTATGCTACTACAAAAAGTTGTTAACT
AAATTAACATTGGGAAATCTTATATTCCATATATTAGCATTTAGTCCAATGTCTTTTAAAGCTTATTTA
ATTAATAAATTTCCAGTGAGCTTATCATGCTGTCTTTACATGGGGTTTTCAATTTGCGATGCTCGATTA
TTCCCTGTACAATATTTAAAAATTTATTGCTTGATACTTTTGACAACAAATTAGGTTTTGTACAATTGAA
CTTAAATAAATGTCATTAATAAATAAATGCAATATGTATTAATATTCATTGTATAAAAAATAGAAGAA
TACAAACATATTTGTTAAATATTTACATATGAAATTTAATATAGCTATTTTTATGGAAATTTTTCATTGA
TATGAAAAATATGATATTGCATATGCATAGTCCCAGTGTAAATCCCATTCAACTTTCATTAAAGCA
TTTACTTTGAATTTCTCCAATGCTTAGAATGTTTTTACCAGGAATGGATGTCGCTAATCATAATAAAT
TCAACCATTATTTTTTTCTGTTTATAAATACATTGTGTTATATGTTCAAATATGAAATGTGTATGCACC
TATTGAAATATGTTTAAATGCATTTATTAACATTTGCAGGACACTTTTACAGGCCCAATTATCCAATAG
TCTAATAATTGTTTAAAGATCTAGAAAAAATAAATCAAGAATAGTGGTATTTTTTTCATGAAGTAATAAAAA
CTCGTTTTGGTGAA
ACGCGTAAGCGGCCGCGGCATCTAGATTGCAAGAAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTTGATTCCACCGCCGCTTCTATGAAAGG
  
```



[View online »](#)

Restriction Sites:	Sgfl-Mlul
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	NM_080921.4
Summary:	The protein encoded by this gene is a member of the protein tyrosine phosphatase (PTP) family. PTPs are known to be signaling molecules that regulate a variety of cellular processes including cell growth, differentiation, mitosis, and oncogenic transformation. This PTP contains an extracellular domain, a single transmembrane segment and two tandem intracytoplasmic catalytic domains, and thus is classified as a receptor type PTP. This PTP has been shown to be an essential regulator of T- and B-cell antigen receptor signaling. It functions through either direct interaction with components of the antigen receptor complexes, or by activating various Src family kinases required for the antigen receptor signaling. This PTP also suppresses JAK kinases, and thus functions as a regulator of cytokine receptor signaling. Alternatively spliced transcripts variants of this gene, which encode distinct isoforms, have been reported. [provided by RefSeq, Jun 2012]
Locus ID:	5788
MW:	51.8