

Product datasheet for SC213400

SLITRK3 (NM_014926) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	SLITRK3 (NM_014926) Human 3' UTR Clone
Symbol:	SLITRK3
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_014926
Insert Size:	1207 bp
Insert Sequence:	>SC213400 3'UTR clone of NM_014926 The sequence shown below is from the reference sequence of NM_014926. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAACGATCGCC
GTCCTGGAGAAGACAACATACAGGTTCTAACAGAGAGAAGAAAATATATTAGTGCTTTTTTTTTTCAA
AAGAAAAGGAAAAATAAAGAAAATATATCCCTTGCTCCCTTTACTTGTCCAGTAACCTCATCTCAC
GATCTTTCTACCTGAACAAAACCTAAACCGCATGATAACTAGAGAATACAGATGTATGCTCTCCCT
CTCAGATGCGATTGGAGGAAGGGCCATACTCAGATCATTAAATCAATGAAAGTGCCTTCGCAGACTTTT
GCCAGCAAATGTTATCATTATTTTTTTTATACTGAAACTTGAGACTTTGACTGTGCCATGTATAAGATAT
ACTGGGGATCATTGTATGGATCCTAATTAAGTAAATCAATGTGTCTTTTTATTTTCACTAACTATTT
TTTTTATAGTTGTAGTTTTGATTTAAAGGGGGGAAACAAGTTGACATTTGTCTTTGTGGCTTTCTTT
CTTATCATCATGGCACAGATTCTGTACATGTATTAACAATGCAGTTTGCTGCATGCCTGGAACTGCAA
GATGGGGAGGGAGGGGGCGGGTCTTGCTCGAATGTTCTCACTCACTATCTTGTCTCTACACACATATC
CATCTGCAGACCATGGTCCCTAAACCTGCAGTTTCTTCTGCTGGTGAGGTACACACACACACACACA
CACACACACACACACACACACACTCCATACCATTTCATTCACTCAATCTGCTTAGTTCGGGTTTG
ATCCATTTTTCTCTCTCCATAGTTCTACCACCCTTCACTGCTAGTGTACAGCTCACAGCATCTCTCT
ACCACCCTGGGAAAAGTACATTCTAGGCTGGATTCTACTGAAGTAGAGGCTGGTGGTTTAAACAGCCG
AAGACACACGCCTAGGGATGAGCACCTCTTTGTGACACTTCATCCTGATGCCAAGATTTTTTGGAGAA
CATCTGCACTTTCTTATATATATATAATATTAATAAAGCAACTGGGTATATACCTAACAGCTTT
GTAGGAAGCACTTTTAAATGTTTTCTCTCTCACACAAAGATTCAAAAGAAATGTGCATATATTTATTC
TGTAATTTCACTGATTATAAATGTAAAGGTAATGTTTAAATCATTGTATAGTGATTATGCCTCTGGTAT
AGCTTTCTTAATAAAGTTTTTGAACATAA
ACGCGTAAGCGGCCGCGCATCTAGATTCTGAAGAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
  
```

Restriction Sites: SgfI-MluI


[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_014926.4</u>
Summary:	This gene encodes a member of the Slitrk family of structurally related transmembrane proteins that are involved in controlling neurite outgrowth. The encoded protein contains two leucine-rich repeat (LRR) domains and a C-terminal domain that is partially similar to Trk neurotrophin receptor protein. Enhanced expression of this gene was found in tissue from several different types of tumors. Alternative splicing results in multiple transcript variants, all encoding the same protein. [provided by RefSeq, Jan 2016]
Locus ID:	22865
MW:	46.9