

Product datasheet for **SC213402**

TMCC2 (NM_014858) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	TMCC2 (NM_014858) Human 3' UTR Clone
Symbol:	TMCC2
Synonyms:	HUCEP11
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_014858
Insert Size:	1249 bp



[View online »](#)

Insert Sequence:	>SC213402 3'UTR clone of NM_014858 The sequence shown below is from the reference sequence of NM_014858. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC CTCCTGGAGCACGTGTTGCTGCCAGCTGAGTGGCCAGCCACCAACCCTGTGCTCTCTGGCCCCAG CTGGCCACACTTCTCCAGGAGGGACCTTGGACTTCTTTGTGTGTCAGTTTGGCCTCCTGCCCAAAT GTCCATTCCAGCAGCTCCTGCCCCCTTCTCTGTACTTGCTTCTGTCTGACACCTTCTCCCTGTTGGCCT GAAGGGAGCTTAGAATGCAGCCCTACCTGGAGATAGTGCGGGCACCTGTGGCCAAGTGGAGCAGAGGTG GACATGGGGTTGGATTGTTTTGATTATTTATAGTTACACAAGGACTTCTCCAGCTGACCCTCAGGATG CCCCAAGTCAGGAAGACCATTAGAATAGGAGGAGAGGGCTCTGCCTCACTTTCTAGGAAAGAGCCC ACCTCGGAGATAGCTACGGTTTTCTCTGGTGGAGATGGTGAGGATGAAGGCTGGAGAGTGAGGGAGGAG GCTCTGCTGGCCGAGAGAACACAGGGATGGGAGGGTCCCTAGCCTTCGGGCACCTCCAGGGCCAGAGA GCAGGCTCAGAGCAGCTAGTGTGGAGCTCAGCATCCCCACCCACCCCTCCTCCCTGTAGAGCTGATTT GAGGCCTCCTTCTGGGGCTGGGCTCTGCAGGCCAGGTGGGTGTGGCCTGTGTTTTCCCTTCTGTTCTTT CTGCCTGTACTGGATCTGTTATTTTTCAGGGAACAGGCCCCAGGGCCCCCTGAGCCTCACCCCTAAGCC CTTAGGCCTCTGAGAGTGCTGTTGGGTTCTATTTATTTATTTATTTGTTCTTTGTTCCCTACCCGTGC CCCCAGTGTCTTCCCTGCTGAGTACCAGGAGAGGTCTGCCCATCCTCTCTCTGAAGCCAGGGCCCTT CCATTCCATTTAGCCTTTGGATCATCTGGCTGGGAGAAGTGGGACCGAGCCACCCAGCCCCACTATCC CCAAGCAGCCCTACAGCCGGGATGGGAGGCACGTGGCCTCTCTTTATCCGTCTATTTATTTTGTAAAGT GTATTCGTGTGGAGGAGTTGTTGCTTTATTTTTTAAGGCTCTGGAGTGTGTGTATGGTTTCTTTTC ACATCCCAGCCTCCCATGGGCACCTCTAAGAAGAGAGGGGATTTCTTGAAAAGGAGAGAGGAATCCCC TAGAGCAGGGAAGCAGTGCCTGCCAGCTGTTGTGCACCTTCTGAGAAATAATATCTCTAAATTTT CAAACCA ACGCGTAAGCGGCCGCGGCATCTAGATTCAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
Restriction Sites:	Sgfl-MluI
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_014858.4</u>
Summary:	May be involved in the regulation of the proteolytic processing of the amyloid precursor protein (APP) possibly also implicating APOE.[UniProtKB/Swiss-Prot Function]
Locus ID:	9911
MW:	46.1