

Product datasheet for **SC213411**

RIL (PDLIM4) (NM_003687) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	RIL (PDLIM4) (NM_003687) Human 3' UTR Clone
Symbol:	RIL
Synonyms:	RIL
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_003687
Insert Size:	1255 bp



[View online »](#)

Insert Sequence:	>SC213411 3'UTR clone of NM_003687 The sequence shown below is from the reference sequence of NM_003687. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site <pre> GACGCCCCGAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTGTAACAATTG GCAGAGCTCAGAATTCAAGCGATCGCTTGGCGCGCC TACCCCAATGCCAAGGTGGAACCTCGTCTGAGCTGGGACCCTGCTCCACGCCTGCTTCTTAAGGTCCCT GCTCGGCCGGTGTAAATATGTTTCACCCTGTCCCTCTAATAAAGCTCCTCTGCTCCACCTTGAACCTGT CACCTGGCCTGCCACCCTCCTGCGCAGGCCATGCATGGCTCCCGAGTACAGTAGTATCTGCTTAGGTGC CAGGCATGTCTAGGCTCTGGGTGCAGTAGTGAGCAGGACGGGTACCATGCTGCCCTGAAGGGGGCCAC AGCCTGGGTGGAAGGCAGACCTGAATCACAACGGGCCAGCTCTAGTAATACGAAGGTGAGGTCTGGGAT GCTGCGTGC GGCGCAGCAGGAGGTATGACCTGGGTGGGTGAGGAAAGTCTTAGCTGAGAACTAA GAGATGAGGGAGGCACAACTCTGCACGGGGAACTGTGTGTCAAAGGTGAGCTGGGGGCGAGAAAGG CCTCTGTGGCCGTCATCGGCACTAGCGGTTGGGGGTGGTGGTGGGTGCTGGCAGCCTCAGGAGCGAT CGCTCCCCTGGAGAGGGTGATTGAGTGTGAGATCACCGGGGCCCGCGGGCCGTCTGCACTTCTGTT CCAGGTTTTCTGGCATTCTGTTCCAGGTTTCTTCCACCCACGAGCACTTTATTCTCCTCCGAGACC CCCTTTTCTTCCCTTCTCTCTCCCCACCAGGGGGCGCGCTTAAATTTCCAGACAGAGACCAGAAGG AAGGTTTAGAAAGAAAAGAACACTTGCCAGGGCAGCTTTGCCCTCCAAGAGGCTCCCTCCTGTCTCTA AGTCAATTCCACCCCTCCTCGAAGTCCCGGCTTCCACCCCTCCTCAGAGTCCCGGCGTTTACCTCGT GGTGTGTTTGCCTTGGGCCTTTCATGCCCGGCTGCACACATCGTCTGTGAATTGTGGCTGTTCACTCC CGAGGATGTGTCTAGACTCCGGGTGCGTGGATCTACCCTCTAGTTTACTTGCTCGGGAGAAGAACT GACTCGTTTTATTAGTGCCTATTTAGCGAGCCAGAGTAACGTACATTTGTGCTGTTTTCAATTTTGT GCTATCGCAAATCACAAAAAACTGTTATCAATTACATTATCATCAGCACGAGACCCATTTCCTTGT GCTCCTGCCAGCTCAGGATATTATGTTTCTTTTTCATTTTGGCAGTCTGTTAAGTAAATGTCATAT GCGTTTGTTTTAA ACGCGTAAGCGGCCGCGCATCTAGATTCAAGAAAAATGACCGACCAAGCGACGCCAACCTGCCATCA CGAGATTCGATTCCACCGCGCCTTCTATGAAAGG </pre>
Restriction Sites:	Ascl-MluI
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	NM_003687.4
Summary:	This gene encodes a protein which may be involved in bone development. Mutations in this gene are associated with susceptibility to osteoporosis. [provided by RefSeq, Nov 2009]
Locus ID:	8572
MW:	45.2