

Product datasheet for **SC213419**

ZNF384 (NM_001039920) Human 3' UTR Clone

Product data:

Product Type: 3' UTR Clones
Product Name: ZNF384 (NM_001039920) Human 3' UTR Clone
Vector: pMirTarget (PS100062)
Symbol: ZNF384
Synonyms: CAGH1; CAGH1A; CIZ; ERDA2; NMP4; NP; TNRC1
ACCN: NM_001039920
Insert Size: 1267 bp
Insert Sequence: >SC213419 3'UTR clone of NM_001039920
 The sequence shown below is from the reference sequence of NM_001039920. The complete sequence of this clone may contain minor differences, such as SNPs.
 Blue=Stop Codon Red=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
ATCCAGGTGGAGCACCTGGCCAGCTCTTAGAGATCCGTGCTGCCACCCACTGGGAAGAGGAAGAAGTAG
TCCTGGTGTCTTCTTTCTCCAACCTTTGGTGGGAAAAGTCCTTTTCTCCTTGACAGGCCTGGCTCCA
TCTCCTGGGCCTCTGTACGGCTTTCTTCCACAGGATACCATCCTTTTTCTGAACTCTTCTCAAAG
GAACATCAGCCCTCTGATTGCAAAGGAATACTGAGCTGATGGTGCATCCAGCAGCCTCCCTCCCAA
GCAAAGCTTCTAAAAGTGGGGTGGTGTCAAGGGAAGGATTTGCTATGACCTCATAGAACCTTGTC
AGTGTGGCCACTTACCCTATCCTTACCCTCCTTATCCTCAAAGTTTGGGCTGATGTAAGACTAGAGGCT
GGCCCTCCAGATAACAGAGAAAAGGGAGCCCAAATGCAACCAGCCTCTGTTCTATTCTTGCCTGCA
AAAGAACAGAGGTTTCTCAAATGCCTCAGTCCCTGAGAGCCATTTCTCCCTACATCGTCTCACTTTG
CTTCTATTGACTGCTGGTGAAGGAGATTTGGGGTAGGGGCTAGACCTCCTTTTATTTGAAGGGGCA
AGGGCTGAGATGTGGTCCCAAGGGGCCAGAAATCCCAAGTTGGTACAGGTGGCTTAGAAGTGTGTG
TTATGGTTTTACGGATTTCTTGAAGCCTCTCCTTCTCTGCCTACAAAGACCTATACTCTCAGTCT
CCCAACCCACCCCAAGGAGCTGTGGAGGCTTTGTGTTATCTGTGAACTCCAAAACAGGGGTGTG
CGGAGAAGGGAGAGTTCAAGGCAAACGCAAGGACTGGACTTAGCTCCCTAGGTGCCACAGTCAGATGCC
GGACACGGATTTATATATAAATATATATAAATATATACTACTCATCAGGCCATCTTTGTT
GTAACCATTTCTGTGTTATAAATGCATTATCTCTGAGAATTTTCATATTTGATGTTTTGTTATTTTT
GTCCTTTTTTCCCTCTCCACCCTGTCTCTAGCCACAGCATTTTTCTTTTTGTCTTTTTTTTTTTT
TTTTTAAATCATGGCAGATTTAGAGGAAAAGAAAATTAATAAATAAATCAGGAAACAGTTGTTATAAA
GTAATTTAAAAATGAAGAAAAAAGAAAAAATTAATGTACAAACCAAGGGGTGTTTTAGAACATTGT
ATAGAAAATAAATTCATGTAAAAGGA
ACGCGTAAAGCGCCGCGCATCTAGATTGAAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCCTTCTATGAAAGG
```

Restriction Sites: SgfI-MluI



| | |
|------------------------|---|
| OTI Disclaimer: | Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs). |
| Components: | The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials. |
| RefSeq: | NM_001039920.3 |
| Summary: | This gene encodes a C2H2-type zinc finger protein, which may function as a transcription factor. This gene also contains long CAG trinucleotide repeats that encode consecutive glutamine residues. The protein appears to bind and regulate the promoters of the extracellular matrix genes MMP1, MMP3, MMP7 and COL1A1. Studies in mouse suggest that nuclear matrix transcription factors (NP/NMP4) may be part of a general mechanical pathway that couples cell construction and function during extracellular matrix remodeling. Alternative splicing results in multiple transcript variants. Recurrent rearrangements of this gene with the Ewing's sarcoma gene, EWSR1 on chromosome 22, or with the TAF15 gene on chromosome 17, or with the TCF3 (E2A) gene on chromosome 19, have been observed in acute leukemia. A related pseudogene has been identified on chromosome 7. [provided by RefSeq, Apr 2011] |
| Locus ID: | 171017 |
| MW: | 48.1 |