

Product datasheet for **SC213431**

AFAP1L2 (NM_001001936) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	AFAP1L2 (NM_001001936) Human 3' UTR Clone
Symbol:	AFAP1L2
Synonyms:	CTB-1144G6.4; KIAA1914; XB130
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_001001936
Insert Size:	1248 bp



[View online »](#)

Insert Sequence:	>SC213431 3'UTR clone of NM_001001936 The sequence shown below is from the reference sequence of NM_001001936. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC AAGGAATGGGAGAAGAAAGGAGCAAGTTAAGAAAAACAAGCTTCATCTAAAGACTCTCATGTCAATGTGGA CCTTGGTGACAATCCTGCTTTGTTAAAGCAAAAACTATGCGAAAGGGTGAGTCTGTTTAGAAGAAAAAG CAAAGACTGAGGTACTGTGAATGGAGAGCTTCAGCTAAGAGGAGGCTCTGTCCCTTTTCAGAGCCAAAG GAAATAATACAACAAAAAGGAGGCTTCTTTGGAGACCTAAGTCTATTGGATGTAAACAAGACGTTGTAT TTAGGGATGTTCTGTGTTTCTTTCTTTTGAAGTTGTCATCAATTGCTTTACTAAGATTTTAAATAG TGAAAACCTCCTGTTTAGACTTTGGTGAAGATGAATCAAGGAAGCAGGGCCCTGTCTTATGGGTCACG TGCTTTTGGTGAGTGAGAAGACCTAAACTCCTGGCCATCATCTCTATCCAATACTTAGCAGTTGGGGA TTAAACCATCCTTGCCTTCAGTTCTCTCAATATTACCAGGCCCAACTCAGTCTTCAGTGATTTTAAAC AGCATTGACATCATCTGTAAAACCATCATCTGTAAAACCATCTATGACATGAGTTTTGAGAAACAATAA TGGGGAAAATATTTGGGACCAAGCTGAAGCACTAATCCCACTAAGTTAAAGACTTCTTTCCAGTCCAAG GCAGGCCCTGAATCAACTGTCTTTAAATAAAATTTTAAAGTGATGCTGTATTATATATAGGAAAAAATGCT TAAATCCTGTCAATTTAGAACAGTGAAAAGTATCTTTGAGATTAAGTGACTCTTACTGTAGGAAAA ATATTACTCTGTGTTTACAGATTCATTGCTGTGGTCAGGCCATTTTAAAGGAAGAGTTATTTAATATA AATAGTCTCTGATTTTAAAGTCTGTTTAAATGTTTCTCTTCCAAGAACAAAGTGGTGATTTTGGT TAGGGTGATCGCCCTCTTAAATTTGGCAGTGCTGTTCTTGTGCTGCCCTGTCTTTCTCTGTATGGC ATTTTTTTTTTTTTTTTTTAAACACAGGTTGAAACATTTTCATCTATTATCTCTGCCTCATTCTGGAGG GTTGTGATCAGTTCTCTAACACTTGTTCCTGAGAACTAAATGTCTTTTTTATTCTTATTTCTCTCTC ATAAACATTTGGTGACCTTTTACCAAGTGGTGAGTTAGTTAGTTTTTAAATAAAATGTTCAATTGTA TTTGAA ACGCGTAAGCGGCCGCGGCATCTAGATTGAAGAAATGACCGACCAAGCGACGCCCAACCTGCCATCA CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
Restriction Sites:	Sgfl-MluI
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	NM_001001936.3
Summary:	May play a role in a signaling cascade by enhancing the kinase activity of SRC. Contributes to SRC-regulated transcription activation.[UniProtKB/Swiss-Prot Function]
Locus ID:	84632
MW:	47.8