

Product datasheet for SC213457

ARMC7 (NM_024585) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	ARMC7 (NM_024585) Human 3' UTR Clone
Symbol:	ARMC7
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_024585
Insert Size:	1264 bp
Insert Sequence:	>SC213457 3'UTR clone of NM_024585 The sequence shown below is from the reference sequence of NM_024585. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAACGATCGCC
CCGAGGAGCGTGGCCCCACGGCAGCGCTTGATCCATGGAGACTGCGAGACCGTGGCACCCCTACTGCTGG
AGACCACAGTCTGATGTGGACGCAGGGAACGGGGAGCACATACTGCCCCATTGGTGCCTTTTCAGCCA
TCTGAAAGCGGGTTCTTTCAGCAGGACAGGCATTTACACTGATGAAACGCCACTGGGAGTGAGGAAGC
CAGACTCCAGAGACACGGAGAAGATCAAACCTGGAGCTGCGTTTCATAGGCTGGCACTCTCAATCCTACAT
CAGGTGCCACCAACCAGACTCAGGCCCTGGTGTAAAGAAGCGGCCAAGTGCCTGGACCCAGAGGCTTT
GCAGGACAGTGTTCTCAGGAGCTGGGCCTGAGGCTTAGGAGAGCTGCCTTCGCTGCAGGAAATCAGGGA
TTATCCCTTAACAGAAGTGTCTGGAGTAGTTTTCAGGTATAGGAATGAGATGCCTCGTGGTGAAAGGAT
CTCACCCTGGGAAGATGTGGTGGCCCTCCAGGGCTCTGGAGGATGGATGCCTCCCCAGGGGCTCTCC
AAGCTGGGCATTTGGGCCTGGTGGATGCCAACCTGGATAACCTGTGGCCAGCATTGACTGTCCACCCA
GCCTTGCTGTTAGGCACCATGACTCCAAGATGAAGATGTGGTCCCTGCCCTTGAGTGACAGCCAGGGA
CTTAATGTGGCCATCGGGCATCAAGCACAAGGCCATGCAGGTGATGATACGTGGAATAGAGGCACCAG
CCCTGGTAACTGCATCTTCTCCCTTGCCACCCATGGCCCCGGCTGAAAGCTTCGGCCCTCCTCTGCT
GTCACTCAATGATGGGGAGCCCTACCCAGAAGTGATCCACGAGGGCATCAGGGACGCAGTGAGTGT
TGCTCAAGGGAGTCAGGAAGAGACGGCAACGTAAGGATGTGGCTCCATGTCCATGGTGGCCCTGGTC
AACATAAGGAGCGTGGGATCCGATGGAAAGGTGGAGCTCAGGGAAATGGGGGTCTTGCCTCTCGTGT
ACCCCTCAAGGCTGACCCCTTAGATGGCCAGGAATGGCAGGTGCTACAAAAATGGTACCCACGTGGG
CATGGAATGGGGCAGATTAGGGGACCACTGGACTCAGAGGGGAGGGAAGGGTCATCAGCACCCGCTC
AGGGAGCTGTCCCTTTATGTTCCCAAATAAAGGGTCCTAGAAGACTAGAAAGCCAAGGTCTTTTATTA
AAGGTCCCGACTGGTTTTCCCA
ACGCGTAAGCGGCCGCGCATCTAGATTGAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
  
```


[View online »](#)

Restriction Sites:	Sgfl-Mlul
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_024585.4</u>
Locus ID:	79637
MW:	46.1