

Product datasheet for **SC213459**

IPP (NM_005897) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	IPP (NM_005897) Human 3' UTR Clone
Symbol:	IPP
Synonyms:	KLHL27
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_005897
Insert Size:	1285 bp



[View online »](#)

Insert Sequence:	>SC213459 3'UTR clone of NM_005897 The sequence shown below is from the reference sequence of NM_005897. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site <pre> GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG TAACAATTGGCAGAGCTCAGAATTCAGCGATCGCC CGTTGTGAAGGGGCGTTGCTGTGCTATGAATATAAAGCATAACAGAGCACAAGGAACATTTTGAAAA TGCAGGACCTGACCTTTTGCAAAGCAAACATAGATTTGGTGATAAGACCTCTGATTCTGTGAGTTTTC CATGTATTTGGTAGATAAATGACTAACGAGGCTGGGCGCAGTGGCCTCATGCCTATAATCCCAGCACTT TGGAAAGGCTGAGGCAGAAGGATCACTTGAGCCCAGGAATTTGAGGTTGTAGTAAGCTACGATCATGCCA CTGCACTCCACTCTGGGCAACAGAATGAGATCCTGTCTCTTAAAGAAAAAAATAAACTCAAGAATGAT TTTTCTAAAAATATGAGTGAGAAATGGATATGGAATATATCCCAATTAATAAGAGCTCACACAGATTG CCAGGTTTACTAATAATAGAAATGCCATTTATGCTGTAGAATTCTCCAGTTTAAGGTAGACAACCTCT ATTTTACAAATTTGTAGAAATTTTACTAGAAGCTTCCTGGCAAATAGTCCCTTTAGCTACTTCTCCA GTGCAGTATAAATTGTAACATAAACTTTATGTGATATAGGAATAGTGCAATGACACAGTTTAAAGAC CTCTATATTCACAATTATATCCAGAGCATAATTTTCCATGGCCATAGTTTGAGCCACAGAGAAAGCCAA TGAAATAGATTCTACTTTTATTTAATGCCAAGAGTTTCTACTAGCCCCAACTTTTATTGGCAGCAG ATATCCTAGGAGGGGTGTTTGGGAAATATCTAGCTATTTTAATAGATAGATACCTTTTGTCTGATGAT AAGTTTTTGCTTGCTCATATTTTGTCTTGATCCTGAGAATGCAAGTAAGGAATTAATAGGAGTTCT TAGTGACCAGATGTTGTAATGGGCTATACAGTCCAGCAGGAGTATGCACCTAGAAATTTCTCCCTTTT TTAGGCTGGGCATGGTGGCTCACACCTGTAATCCCAGCACTTTGAGAGGCCGAGGTGAGACGGGAAGAT GGCTTGAACCTCAGGAGTTCGAGACTAGCCTGGGCAACGTGATGAGACTCCGTCTCTACCAAAAAATATA TAAAAATTAGCCGGGCGTTGTGGCTTGACCTGTGGTCCCATCTGCTTGGGAGGCTGAGGTTGCAGGAT CACTTGAGCCCAGGAGGTCAAGGCTGCAGTGAGCCAAGATTGTGCCACTGCACACCGGCCTGGGTGACA GCATGAGAACCTGTCTCAAAAAAAAAAAAAAAAAAGAAAAA ACGCGTAAGCGGCCGCGCATCTAGATTGAAGAAATGACCGACCAAGCGACGCCAACCTGCCATCA CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG </pre>
Restriction Sites:	Sgfl-MluI
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	NM_005897.3
Summary:	The protein encoded by this gene is a member of the kelch family of proteins, which is characterized by a 50 amino acid repeat which interacts with actin. Transcript variants have been described but their full-length nature has not been determined. [provided by RefSeq, Jul 2008]
Locus ID:	3652
MW:	48.6