

Product datasheet for **SC213467**

SAM68 (KHDRBS1) (NM_006559) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	SAM68 (KHDRBS1) (NM_006559) Human 3' UTR Clone
Symbol:	SAM68
Synonyms:	p62; p68; Sam68
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_006559
Insert Size:	1283 bp



[View online »](#)

Insert Sequence:	>SC213467 3'UTR clone of NM_006559 The sequence shown below is from the reference sequence of NM_006559. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site <pre> GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC TACAGAGAGCACCCATATGGACGTTATTAACAAACAAACATGAGGGGAAAATATCAGTTATGAGCAAAGT TGTTACTGATTTCTTGATCTCCAGGATTCTGTTGCTTTACCCACAACAGACAAGTAATTGTCTAAG TGTTTTCTTCGTGGTCCCTTCTTCTCCACCTTATCCATTCTTAACCTCTGCATTCTGGCTTCTGT ATGTAGTATTTTAAATGAGTTAAATAGATTTAGGAATATTGAATTAATTTTTTAAAGTGTAGATGC TTTTTCTTTGTGTTTAAATATAAACAGAGGTACCTTTTATAATAAAAAAAGAAAGTTGAGTAAAA AAAAAAACACACAAACCTGTTAGTTTCAAAAATGACATTGCTTGCTTAAAGTTCTGAAGTAAAGGCT TGTTAAGTTTCTTAGTTTGTGTTGAGGCATCCCGTAAAGTTGAGTTGCAGAATCCCAAAGTGGC TACATTTCAAAATTCAGGGCTGTTAAGATTTAAATCACAAACATTAACGGCAGTAGGCACCACCATG TAAAGTGAGCTCAGACGCTCTTAAAAATGTTTCCTTTATAAAGCACATGGCGGTTGAATCTTAAGG TTAAATTTTAATATGAAAGATCCTCATGAATTAATAGTTGATGCAATTTTAAACGTTAATTGATATAA AAAAAAACAACAAATTAGGCTTGTAAGTACTGCTTTTCTTACGTGGGTTTGAATCTAGCCCC AGACATACTGTGTGAGAGATACTTAGAGGGAGGAGTAGGTTTTGAAGAGGTTGATGGTGGTGGGAG GGAAGGCTCCTGAATTGAGTTTGTGAGAGCTTTTAGCCATGAAGAATCTTTCAGTCATAGTACTA ATAATTAATTTTTCAGTATTTAAAGACAAAGTATTTTGCCATTGAGATTCTGCACTCCATGAAAA GTTCACTTGGACGCTGGGGCCAAAGCTGTTGATTTTCTTAAGTTGACGGTTGTCAATATATCGAACTG TTCCCAAGTTAGTCAAGTATGTCTCAACACTAGCATGATATAAAGGGACACTGCAGCTGAATGAAAA AGGAATCAAAATCCACTTTGTACATAAGTTAAAGTCCTAATTGGATTTGTACCGTCTCCCATTTTGT CTCGAAGATTAAATGCTACATGTGTAAGTCTGCCTAAATAGGTAGCTTAACTTATGTCAAAATGTCT GCAGCAGTTTGTCAATAAAGTTTGTCTTTTAAATCAAA ACGCGTAAGCGGCCGCGGCATCTAGATTGAAGAAATGACCGACCAAGCGACGCCAACCTGCCATCA CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG </pre>
Restriction Sites:	Sgfl-MluI
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	NM_006559.3
Summary:	This gene encodes a member of the K homology domain-containing, RNA-binding, signal transduction-associated protein family. The encoded protein appears to have many functions and may be involved in a variety of cellular processes, including alternative splicing, cell cycle regulation, RNA 3'-end formation, tumorigenesis, and regulation of human immunodeficiency virus gene expression. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Dec 2012]
Locus ID:	10657
MW:	48.6