

Product datasheet for **SC213487**

FAM61B (LSM14B) (NM_144703) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	FAM61B (LSM14B) (NM_144703) Human 3' UTR Clone
Symbol:	FAM61B
Synonyms:	bA11M20.3; C20orf40; FAM61B; FT005; LSM13; RAP55B
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_144703
Insert Size:	1256 bp



[View online »](#)

Insert Sequence:	>SC213487 3'UTR clone of NM_144703 The sequence shown below is from the reference sequence of NM_144703. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site <pre> GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC TCCCACAGGGCCGGGACTGGCAGGGTGAGGGTGCAGCCAAAGGCTCCTACTGAAGTGGCGCATAACT GACGCTGTGTGTGTGACGACGCGAGGAAACGCTGCACCTACAGGGAGAGGTGGTCACTTTGTTTACGG AGTTTGGAAGAGACCCATACTGCTACTTGTGTTTTGGACTTAACTGAACTTGGACATGGTCTGAGTTAG AACCACTGTGTTTTGGGAAGTATTCATGGGTAACCTCTTTGAGGTCTCTTTATCTGTGTTTCCTTTTAA GTTGCGCATAGCCTAATTCTAAGGTTTTGGTATTTTGCAGGAGGTTTCTATAGTGAAGCTGAATCCT TACTTTGTGACTTTTTTTTTTTTTTTTAAATGACAAGCTTTGACTTTTAAAGTGAACCAAAATCTGTTG GCAGAGGTGGCAGCCAAGTACATCTCTGAACCCAGCTGGCCCTGGTGTGTTGGCCTGGCAGCCAC TGCCAAGGGTGGGGTCTCAGGAGTCAGGCAGGGCCAGCACAGGGTGGCGTGGGGGGCAGGGGTGGGTGG GTGGAGGGCAGGAAGGGTTTTCCCATGGATCATGTTGTATAAGTGAACCAGACCACCTGATGGCAT CCACAGTGATGTCAAGGTTGGGGCTGGCCAGGGGTGGGTGGACTAGAAGCATTTGGGAGTAGTGGCCAG GGGCCCTGGACGCTAGCCACGGAGCTGTGCACAGAGCCTGGTGTCCACAAGCTTCCAGGTTGGGGTTG GAGCCTGGGATGAGCCCCGGCAGCGCCTTGGCCCTTCTGTGGTCCCTGCCAGCCTCTGACCTGGGCCGG TCAGTCATTGCTGGACTCTGGCCACACACTGGCGTTCTCATCCACTTGGAAACAAGCCAGTCTTTCTG CAAGGTCAGTTGACCAAGAGCATATTTCCCTCTGTTGTACATCGTTGTTTTGTGTTTGTGTTGTAACA GTGGGTGGAGGGAGGGTGGGGTCTACATTTGTTGCATGAGTCGATGGGTGAGAACTTTAGTATACGCAT GCGTCCTCTGAGTGACAGGGCATTGTCGAAAATAAGCACCTTGGTAACTAAACCCCTCTAATAGCTA TAAAGGCTTTAGTTCTGTATTGATTAAGTTACTGTAAAAGCTTGGGTTTATTTTTGTAGGACTTAATGG CTAAGAATTAGAACATAGCAAGGGGCTCCTCTGTTGGAGTAATGTAAATTGTAATTATAATAAACAT GCAACCTTTAAAA ACGCGTAAGCGGCCGCGGCATCTAGATTCGAAGAAATGACCGACCAAGCGACGCCAACCTGCCATCA CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG </pre>
Restriction Sites:	Sgfl-MluI
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	NM_144703.3
Summary:	May play a role in control of mRNA translation.[UniProtKB/Swiss-Prot Function]
Locus ID:	149986
MW:	45.9