

Product datasheet for **SC213507**

PDE2A (NM_001143839) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	PDE2A (NM_001143839) Human 3' UTR Clone
Vector:	pMirTarget (PS100062)
Symbol:	PDE2A
Synonyms:	CGS-PDE; cGSPDE; PDE2A1; PED2A4
ACCN:	NM_001143839
Insert Size:	1273 bp
Insert Sequence:	>SC213507 3'UTR clone of NM_001143839

The sequence shown below is from the reference sequence of NM_001143839. The complete sequence of this clone may contain minor differences, such as SNPs.

Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAACGGATCGCC
AATGGCTGCTGCAGCCTTGATGCTGAGTGATCCCCCTCCAGGGACACTTCCCTGCCAGGCCACCTCCCA
CAGCCCTCCACTGGTCTGGCCAGATGCACTGGGAACAGAGCCACGGGTCCTGGTCTAGACCAGGACT
TCCTGTGTGACCCTGGACAAGTACTACCTTCTGGCCTCAGCTTCTCGTCTGATAATGGAAGCAAG
ACTTCCAACCTCACGGAGACTTTGTAATTTGTTTTCTGAGAGCACAGGGTGACCAATGAGCAGTGGGC
CCTACTCTGCACCTCTGACCACACCTTGGCAAGTCTTTCCAAGCCATTCTTTGTCTGAGCAGCTTGAT
GGTTTCTCCTTGCCCCATTTCTGCCCCACCAGATCTTTGCTCCTTTCCCTTTGAGGACTCCCACCTTT
GGGGTCTCCAGGATCCTCATGGAAGGGGAAGGTGAGACATCTGAGTGAGCAGAGTGTGGCATCTGGAA
ACAGTCTTAGTTCTGTGGGAGGACTAGAAACAGCCGCGGGGCGAAGGCCCTGAGGACCACTACTAT
ACTGATGGTGGGATTGGGACCTGGGGATACAGGGGCCCCAGGAAGAAGTGCCAGAGGGGAGCTCAG
TGCTCTGCAGAGAGGGGCCCTGGGGAGAAGCAGGATGGGATTGATGGGCAGGAGGATCCCCGACTGG
GAGACAGGCCAGGTATGAATGAGCCAGCCATGCTTCTCCTGCCTGTGTGACGCTGGGCGAGTCTCTT
CCCTGTCTGGGCCAAACAGGGAGCGGGTAAGACAATCCATGCTCTAAGATCCATTTTAGATCAATGTC
TAAATAGCTCTATCGCTCTGCGGAGTCCCAGCAGAGGCTATGGAATGTTTCTGCAACCCTAAGGCACA
GAGAGCCCAACCTGAGTGTCTCAGAGGCCCTGAGTGTTCCTTGGCCTGAGCCCCTTACCCATTC
CTGCAGCCAGTGAGAGACCTGGCCTCAGCCCTGGCAGGGCTCTCTCTTCAAGGCCATATCCACCTGTGC
CCTGGGGCTTGGGAGACCCATAGGGCCGGGACTCTTGGGTGAGCCCGGCCACTGGCTTCTCTTTTT
CTCTGTTTCATTCTGTGTGCGTTGTGGGTGGGGAGGGGGTCCACCTGCCTTACCTTTCTGAGTTGCC
TTTAGAGAGATGCGTTTTTCTAGGACTGTGCAACTGTGATATGTTCCCGTGGGCTGACCGCTTTG
TACATGAGAATAAATCTATTTCTTTCTACCA
ACGCGTAAGCGGCCGCGCATCTAGATTGAAAGAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCCTTCTATGAAAGG
  
```

Restriction Sites: SgfI-MluI



[View online >](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	NM_001143839.4
Summary:	Cyclic nucleotide phosphodiesterase with a dual-specificity for the second messengers cAMP and cGMP, which are key regulators of many important physiological processes. Plays an important role in growth and invasion of malignant melanoma cells (e.g. pseudomyxoma peritonei (PMP) cell line) (PubMed:24705027).[UniProtKB/Swiss-Prot Function]
Locus ID:	5138
MW:	45.3