

Product datasheet for **SC213515**

MSK1 (RPS6KA5) (NM_004755) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	MSK1 (RPS6KA5) (NM_004755) Human 3' UTR Clone
Symbol:	MSK1
Synonyms:	MSK1; MSPK1; RLPK
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_004755
Insert Size:	2000 bp



[View online »](#)

Insert Sequence: >SC213515 3'UTR clone of NM_004755
 The sequence shown below is from the reference sequence of NM_004755. The complete sequence of this clone may contain minor differences, such as SNPs.
 Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
CTCTTCCAGTTCTCGGACTCAGTAGCTTAGGCATGGTAGGAGTGTATCAGTGATCCATTGCACCTTTAT
TCCTCAGCATATGCCTGAGCGATCTTTTATGCTTTTAAAAATGTTTCCCGTTGGTCTCATTGGAATC
TGCTCTCAATGATTTTTTTCAGGAAAACCTGTTGGTTATCCTCATTCAAAGCACTGGACAGAGAAT
GTTACTGTGAATAGAGCACATATACTCTTTTTAGCAACCTAGCATGATGCCAACAGACTATTCTTGA
AAGAGCAAAGGTTCTGTAAATTAATTAGGGCTAGATTTGAGCTGCTTGAAGTACAGGTTTTCCAG
ATGTCTGCCAACAAAGAAATGACTCATACTGTGATGATACCTTTTGCCTTGTGGACAATGTGGGT
TTTTGAAATTTGCACCCTTCAAACAATGATTTATCAGAGAAAGGGGTCTGTTTTCAAAAAGATTCTGT
AATGAATTTTATGTGTGGCATATACTTATTTCTTGAGAGAAGATTTAACTTATTGTTTTATTTATG
GTTACATATGATGATAACCTGCTATTATTAACCTTTTTCTAAAAAGTGAAAAAAGATATA
AGAACTCAAGGTCCCATACTCTGTATTCGGATCCATCTGAGATGCATGTAAGCTATGTGTATGTTTT
TAATTTTGCACGTCTTTTCTGGCAATTTGTTTTAATGGTTATTGCAGAATATTAAGGTACATGTCTC
TCTGTTTTAAGTAATATTGCACTTTATAAAAAAGTATGAATAAAGCAAATTTTATAAAGTGCACCTG
TTTAAAGCATTGCACTGTATTTTGGCATTATTTTCTACTTTAAATTTGCTCCTCACATTCC
TCTTCTACTTTGTATGCAACAAGTAGAATGGGGCATGTTGTGTCATGTAGTCAGCCACTTATGCACCAA
TGTGAGGAAAACCTAAAGGGAAATTAACCTAAACACTGTGCTTCATATTTGTACACTGTGTTGTACTAC
AGTGAGGAATTTCTCTGTAGTCATATAATTATGTACATAATATTTTGAATCATACCTATGACTTGTT
TGGAAATTTTTCTGTTAAATTTTAAATCCAGAAAGCATATTTTATAAACTTATGCAGAGCACTTTTATT
GCTCAAAAAGTTCTGAATTCATACAGAAAACAAGTACTATGTGATGAAAACATTTTCAAAAAGATTGCT
GCATTTAAAAATACAATTAATTCGTTCCCTATGCAAAAAAAAAAAAAAAAAAAAAAGAGTGATAGGGT
TTGTGTGTATATTTCTTTTAAAGCCATCATGTGAAATGTCAGGACTGAGAAAAGTGATTTTTTTGCT
TTGTTTACAGGAATCATGCTTAAAGAGATAATGCCCAACATGTTATTTTATAGTTTTTGTAAATGATA
ACCTTTTGAAGCATTAGTTTGGAAATTTGGACATGATATTTATTTTCTTTCTTGAGGTTATAGCAGCA
TTAATTTCAACCAGTTATGAATACTAAAAATATTGCACTCTATGTAATAGTAGTATTTTACTAAA
TCAATGTATATTAATAGCACAGGCAAAAAATGGCAAAATATAAAATATTTTCAAAATATTATATTAT
CCTGTTACATTTTGTCCACAGTGGTTATTGGGAATAATTTTTACACTTTTGTGCGCCAAATGGCT
ACATAATATTGAGTTAATGAAAGTTTTCTACTCCAGCTTTCAGTTATTCCTAGTAATTATGGGAAATT
TTGTGTGAATTATTACCAAATATTGCTAATTCACAAAATGTGATTCCATCAGGCACAGATCTCTAGTA
ATCAAGAAAATGACCATATTTATAATTATATTGGACATATTGGGAACATATTTGGCTCTTCAAGCCATT
TTTTAACTGTGCATCTGGTAGAAAATACTAGTTATAAAGTACAACTAAAGGGACTATGTTGCTAA
ACGCGT AAGCGGCCGCGGCATCTAGATTCGAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
  
```

Restriction Sites: SgfI-MluI

OTI Disclaimer: Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).

Components: The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.

RefSeq: [NM_004755.4](#)

Summary:

Serine/threonine-protein kinase that is required for the mitogen or stress-induced phosphorylation of the transcription factors CREB1 and ATF1 and for the regulation of the transcription factors RELA, STAT3 and ETV1/ER81, and that contributes to gene activation by histone phosphorylation and functions in the regulation of inflammatory genes (PubMed:11909979, PubMed:12569367, PubMed:12763138, PubMed:9687510, PubMed:18511904, PubMed:9873047). Phosphorylates CREB1 and ATF1 in response to mitogenic or stress stimuli such as UV-C irradiation, epidermal growth factor (EGF) and anisomycin (PubMed:11909979, PubMed:9873047). Plays an essential role in the control of RELA transcriptional activity in response to TNF and upon glucocorticoid, associates in the cytoplasm with the glucocorticoid receptor NR3C1 and contributes to RELA inhibition and repression of inflammatory gene expression (PubMed:12628924, PubMed:18511904). In skeletal myoblasts is required for phosphorylation of RELA at 'Ser-276' during oxidative stress (PubMed:12628924). In erythropoietin-stimulated cells, is necessary for the 'Ser-727' phosphorylation of STAT3 and regulation of its transcriptional potential (PubMed:12763138). Phosphorylates ETV1/ER81 at 'Ser-191' and 'Ser-216', and thereby regulates its ability to stimulate transcription, which may be important during development and breast tumor formation (PubMed:12569367). Directly represses transcription via phosphorylation of 'Ser-1' of histone H2A (PubMed:15010469). Phosphorylates 'Ser-10' of histone H3 in response to mitogenics, stress stimuli and EGF, which results in the transcriptional activation of several immediate early genes, including proto-oncogenes c-fos/FOS and c-jun/JUN (PubMed:12773393). May also phosphorylate 'Ser-28' of histone H3 (PubMed:12773393). Mediates the mitogen- and stress-induced phosphorylation of high mobility group protein 1 (HMG1/HMG14) (PubMed:12773393). In lipopolysaccharide-stimulated primary macrophages, acts downstream of the Toll-like receptor TLR4 to limit the production of pro-inflammatory cytokines (By similarity). Functions probably by inducing transcription of the MAP kinase phosphatase DUSP1 and the anti-inflammatory cytokine interleukin 10 (IL10), via CREB1 and ATF1 transcription factors (By similarity). Plays a role in neuronal cell death by mediating the downstream effects of excitotoxic injury (By similarity). Phosphorylates TRIM7 at 'Ser-107' in response to growth factor signaling via the MEK/ERK pathway, thereby stimulating its ubiquitin ligase activity (PubMed:25851810).[UniProtKB/Swiss-Prot Function]

Locus ID: 9252

MW: 77.5