

Product datasheet for **SC213528**

Integrin beta 1 (ITGB1) (NM_133376) Human 3' UTR Clone

Product data:

Product Type: 3' UTR Clones
Product Name: Integrin beta 1 (ITGB1) (NM_133376) Human 3' UTR Clone
Vector: pMirTarget (PS100062)
Symbol: ITGB1
Synonyms: CD29; FNRB; GPIIA; MDF2; MSK12; VLA-BETA; VLAB
ACCN: NM_133376
Insert Size: 1282 bp
Insert Sequence: >SC213528 3'UTR clone of NM_133376
 The sequence shown below is from the reference sequence of NM_133376. The complete sequence of this clone may contain minor differences, such as SNPs.
 Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
GTGGTCAATCCGAAGTATGAGGGAAAATGAGTACTGCCGTGCAATCCCACAACACTGAATGCAAAGT
AGCAATTTCCATAGTCACAGTTAGGTAGCTTTAGGGCAATATTGCCATGGTTTTACTCATGTGCAGGTT
TTGAAAAATGTACAATATGTATAATTTTTAAATGTTTTATTATTTTGAAAATAATGTTGTAATTCATGC
CAGGGACTGACAAAAGACTTGAGACAGGATGGTACTCTTGTCAGCTAAGGTCACATTGTGCCTTTTTG
ACCTTTTCTTCTGGACTATTGAAATCAAGCTTATTGGATTAAGTGATATTTCTATAGCGATTGAAAGG
GCAATAGTTAAAGTAATGAGCATGATGAGAGTTTCTGTTAATCATGTATTAATAACTGATTTTTAGCTTT
ACAAATATGTCAGTTTGCAGTTATGCAGAATCCAAAGTAAATGTCCTGCTAGCTAGTTAAGGATTGTTT
TAAATCTGTTATTTTGCCTATTGCTGTTAGACATGACTGATGACATATCTGAAAGACAAGTATGTTGA
GAGTTGCTGGTGTAAAAACGTTTGAATAGTTGATCTACAAAGGCCATGGGAAAAATTCAGAGAGTTA
GGAAGGAAAAACCAATAGCTTTAAACCTGTGTGCCATTTAAGAGTTACTTAATGTTTGGTAACTTTT
ATGCCTTCACTTTACAAATTCAGCCTTAGATAAAAGAACCAGCAATTTCTGCTAAAAAGTCCTTGA
TTTAGCACTATTTACATACAGGCCATACTTTACAAAGTATTTGCTGAATGGGGACCTTTTGGATTGAAT
TTATTTTATTATTTTATTTTGTAAATGTCTGGTGCTTTCTGTACCTCTTCTAATCTTTAATGTAT
TTGTTTGAATTTGGGGTAAGACTTTTTTATGAGTACTTTTTCTTTGAAGTTTTCAGCGTCAATTTG
CCTTTTAAATGAACATGTGAAGTTATACTGTGGCTATGCAACAGCTCTCACCTACGCGAGTCTTACTTT
GAGTTAGTGCCATAACAGACCACTGTATGTTTACTTCTCACCATTGAGTTGCCATCTTGTTTCACAC
TAGTCACATCTTGTTTTAAGTGCCTTTAGTTTTAACAGTTCACCTTTTACAGTGCTATTTACTGAAGT
TATTTATTAATATGCCTAAAATACTTAAATCGGATGTCTTGACTCTGATGTATTTATCAGGTTGTGT
GCATGAAATTTTATAGATTAAGAAGTTGAGGAAAAGCA
ACGCGTAAGCGGCCGCGCATCTAGATTGAAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
  
```

Restriction Sites: Sgfl-Mlul



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	NM_133376.3
Summary:	Integrins are heterodimeric proteins made up of alpha and beta subunits. At least 18 alpha and 8 beta subunits have been described in mammals. Integrin family members are membrane receptors involved in cell adhesion and recognition in a variety of processes including embryogenesis, hemostasis, tissue repair, immune response and metastatic diffusion of tumor cells. This gene encodes a beta subunit. Multiple alternatively spliced transcript variants which encode different protein isoforms have been found for this gene. [provided by RefSeq, Jul 2008]
Locus ID:	3688
MW:	49.8