

Product datasheet for **SC213533**

UQCRQ (NM_014402) Human 3' UTR Clone

Product data:

Product Type: 3' UTR Clones
Product Name: UQCRQ (NM_014402) Human 3' UTR Clone
Vector: pMirTarget (PS100062)
Symbol: UQCRQ
Synonyms: MC3DN4; QCR8; QP-C; QPC; UQCR7
ACCN: NM_014402
Insert Size: 1295 bp
Insert Sequence: >SC213533 3'UTR clone of NM_014402
 The sequence shown below is from the reference sequence of NM_014402. The complete sequence of this clone may contain minor differences, such as SNPs.
 Blue=Stop Codon Red=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAACGCATCGCC
AATCCAGCTGCCTATGAAAATGACAAATGAGCAACGCATCCGGATGACGGTTCCTGTCTCTGAAAGAC
CTTTCTCTGGAAGAGGAGTCTGCATTGTAGTGTCTCAAAGACACAATAAACTTCCTATGGTCTGCACTG
TTGTGATATTACATTTTTGTGAGTAGAATCCTGTGTGACCACTAATATTCAAGTTTCAAGCGCCCT
CCTCAGCATGAGTCTTGAAGTTTGCCTGGTCTCGGGTTTTGGTAAAGAGCAGGGCTCCCCCTTTGTT
CCATTATCTACAAGAAAGTAGCTGCAGGTACTCCTGCATCTAACCATTTGTGGTGACATGGGCTATGG
CGAGTGTGGTTCTTGGCCCTCTGCTCCCATAGCAACAATAAAGTCTACTAACTATATTGTAATTAGTA
TGTGTCAGGGACGCTTTATTTTTATTTTTGAGATGGAGTCTTGCTCTCACCCAGGCTGGAGTGC
AGTGGTGTGATCTCGGCTACTGCCATCTGTGCCTCCTGGGTTCAAACGATTCTCCTGTCTCAGCCTCC
GAGTAGCTGGGATTACAGGTGCAGGACACCACACCTGGCTAATTTTTGTATTTTTAGTAGAGACAGGT
TTCACCATGTTGGCCAGGCTGGTCTGGAACTCCTGACCTCAGATGATCTGTGCGCCTCAGCCTCCAAA
GTGCTGGGATTACAGGCGTAAGCCATGGTCTCTGGGCCTCAAAACTCTTCTTAGCTACACAAATGAA
CTGTAAACCTCACAAGCACCTTAAACGTATCCTCATTGGCAGATGAGGAAATGGGCACAGAGAGATT
AAGTAATTTGCCAGAAATTATACAGCAGTGAATGGTGGAGTCTTGCACTAAAGTGATGCTCTTAAAACT
ATGCTCTATTGCAGGGGTGTCCAATCTTTAGCTTCCCTGGGCCACACATAAAATACGCTAACATTACA
ATAGCTGATAAGCTTTAAAAAATGGCAAAAAAACCTCAGCGTTTTAGAAAGTTACAAATTTTTGTT
GGGTCGGATTCAAAGCCATCCTGGGCTGCTTGTGCCAGGACTGCGGTTGGACAAGCTTGCTCTACTG
CCTTTCCTCCTAAGGCAGGCCTTAGCCTGAGGGGATATCAGCTATAAAACGGATGAGGTGATCTCTTAG
CCTCCACAGCAGTTTTCTGTTTCTGTGCCAGGAGGCCAACATGTTTTCAACCATATGTGGTTCATAAA
GCCTCAAAGATTTACTCTGGCCCTTTATGAATAAAGTTTACCCCTGATCTAAA
ACGCGTAAGCGGCCGCGCATCTAGATTGAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
```

Restriction Sites: SgfI-MluI



[View online >](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	NM_014402.5
Summary:	This gene encodes a ubiquinone-binding protein of low molecular mass. This protein is a small core-associated protein and a subunit of ubiquinol-cytochrome c reductase complex III, which is part of the mitochondrial respiratory chain. [provided by RefSeq, Jul 2008]
Locus ID:	27089
MW:	47.7