

Product datasheet for **SC213536**

FNDC9 (NM_001001343) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	FNDC9 (NM_001001343) Human 3' UTR Clone
Symbol:	FNDC9
Synonyms:	C5orf40
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_001001343
Insert Size:	1294 bp



[View online »](#)

Insert Sequence:	>SC213536 3'UTR clone of NM_001001343 The sequence shown below is from the reference sequence of NM_001001343. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site <pre> GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC CCCGCTATACTGCCTCATTGTGGGGAATGAGAGTGGGGAGGAGAGCGCCCGCATCATGTAGCCTAAAGC TTTCCACACAGTAGGTCTGTTGTACAAATGTGTCTATAGACTACCCATTTCTCTCCCATCAAACGTAC TGCTATTGTAGGTCACCTGGGTTGGATGAATGCCCATGACAAAGCTTCTCAAGCTGGAAAATGTATCC CCCAGGGTATGCCAAGGCACAAACGAGGCATCTCAGCAGCAGCGTCCACTCGGAGGACACAGTTAAGAA TTGAAAGCATTGTTTTAACATTAACAAACATGGATATGTTATAAAGTAGAAAGCAAAATTTGCATG ACAAATAAACACAGCACTTGAAGGAAACACTTGAGTTACAAGTAGGCTCTCAGAGGTACACTTTTGGG CAAAGGTAGATGTCCAGTTTCTCTGCCCTTGAGGGGATTATTTGAAATAATTTGGGATTATTTGGAA ATAATTTGTGAGAACCCTTGCTCTGGAGTGTTCCTCAACTTTTAGGCAATCACATACCACCGCTCTCTA TTTTTTAAACCATGGATTATCTTCACTATTATTAATAACTTTCTTTAAATTGACTCATTTTTTAA GCGTAACTTATTTTCAAGACAGCCTTATATTACTCCATAAGTGAAAAACCAGCACCTCTCTGCTGCA AACAGAAGGGAAGAGACATAAGAATAAACATCATGAAAAACAAGATAATATTAAATAATCTGACTTAGCT TCTATTGCCTGCCAGTGATTCTGAGCTTAATGCCTGCTTGCTTGTCTTTGATAAAAAGGGAGATCCAGT GTTGGAGAGGTAATAAAGACATATTAGCACCAATTTGTCTTTCTCTTGATATAGCATAAGGTTTGA AAGAGAATCGAAAAGGGAAGTGCTGTTTTTACTTTGCGAGTCAGTGTTAGTTAAATGCCATGTCTGTGA ACCACCTCAAATGCTGTCTAGTCAGCTTCCCTTGATTTGAGAAAATTGTCTTGACCAGACTAGGCAGGCT GTTCTGGAGAAGACGCCTGTACTCTCCAAGCATCATCCAAGACCTTCCCAGTACCCTCTTATATCTATA GAGCACAAAATCCAGAAATCACAGGACGACACTGCAAGTGAATCAACTAAGAAACAGCAGGAGCTGAGAA GCCAAGACGAGAAGCCCCACATCCCTATTCCCTTGCTTACCTCATGCATTCCCTGCTCGGCACCCAGAC TTTTGCCCCATTCTGCTACTTGTGAACAAATAAAGATTCATATACTCACA ACGCGTAAGCGGCCGCGGCATCTAGATTGAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG </pre>
Restriction Sites:	Sgfl-Mlul
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_001001343.4</u>
Locus ID:	408263
MW:	49.4