

Product datasheet for **SC213553**

AIF (AIFM1) (NM_001130847) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	AIF (AIFM1) (NM_001130847) Human 3' UTR Clone
Symbol:	AIF
Synonyms:	AIF; AUNX1; CMT2D; CMTX4; COWCK; COXPD6; DFNX5; NADMR; NAMSD; PDCD8; SEMDHL
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_001130847
Insert Size:	1281 bp



[View online »](#)

Insert Sequence: >SC213553 3' UTR clone of NM_001130847
The sequence shown below is from the reference sequence of NM_001130847. The complete sequence of this clone may contain minor differences, such as SNPs. **Red**=Cloning site
Blue=Stop Codon

CAATTGGCAGAGCTCAGAATTCAA**GCGATCGC**

CTGTGCTCTTGGCAGAAAGGATATA**TGA**AATCAGGGCATCTTTTCAGGAGGATGCTTCCTATTCAGAGGT
TGCAAAGCACTTATGTAAGCCTTAGGTGACTGCATGGGTGGAAGTGCACAGAGTTCACCATCTGAAGCTT
AATCTGGGGCCAGCAGTTCTGATACCTGTTGTATTTGCGTTACAAATGAACTTGAGCACCAGGCCCTCA
TGCCCACTGTCCTCTCGAGCCTTGGGCACAGAAGTGATTCAACTCTCCCGAGAAAGGAAATATGGGAA
AGATCCTCCCGAATACCTCAGCAACTGGACCATGGAAAAAGTGCAGACGAGGGGGTTAAGGTGATGCC
CAATGCTATTGTCAATCCGTTGGAGTCAGCAGTGGCAAGTTACTTATCAAGCTGAAAGACGGCAGGAAG
GTAGAAACTGACCACATAGTGGCAGCTGTGGGCTGGAGCCCAATGTTGAGTTGGCCAAGACTGGTGGCC
TGGAAATAGACTCAGATTTTGGTGGCTTCCGGTAAATGCAGAGCTACAAGCACGCTCTAACATCTGGGT
GGCAGGAGATGCTGCATGCTTCTACGATATAAAGTTGGGAAGGAGGGGGTAGAGCACCATGATCACGCT
GTTGTGAGTGGAGATTGGCTGGAGAAAATATGACTGGAGCTGCTAAGCCGTAAGTGGCATCAGTCAATGT
TCTGGAGTGATTTGGGCCCGATGTTGGCTATGAAGCTATTGGTCTTGTGGACAGTAGTTTGCCACAGT
TGGTGTGTTTTGCAAAAGCAACTGCACAAGACAACCCCAATCTGCCACAGAGCAGTCAGGAAGTGGTATC
CGATCAGAGAGTGAGACAGAGTCCGAGGCCCTCAGAAATTAATCTCCAGCACCCCGGCAGTCCAC
AGGCTCCCGTCCAGGGGGAGGACTACGGCAAAGGTGTCATCTTCTACCTCAGGGACAAAGTGGTCGTGGG
GATTGTGCTATGGAACATCTTTAACCGAATGCCAATAGCAAGGAAGATCATTAAAGACGGTGAGCAGCAT
GAAGATCTCAATGAAGTAGCCAAACTATTCAACATTCATGAAGACTGAAGCCCCACAGTGGAAATTGGCAA
ACCCACTGCAGCCCTGAGAGGAGTGAATGGGTAAAGGAGCATTTTTTTTATTTCAGCAGACTTTCTCTG
TGTATGAGTGTGAATGATCAAGTCTTTGTGAATATTTTCAACTATGTAGGTAATTTCTAATGTTTCA
TAGTGAATAAATTCTGATTC

ACGCGTAAGCGGCCGCGGCATCTAGATTCAAGAAAATGACCG

Restriction Sites: Sgfl-Mlul

OTI Disclaimer: Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).

Components: The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.

RefSeq: [NM_001130847.1](#)

Summary:

This gene encodes a flavoprotein essential for nuclear disassembly in apoptotic cells, and it is found in the mitochondrial intermembrane space in healthy cells. Induction of apoptosis results in the translocation of this protein to the nucleus where it affects chromosome condensation and fragmentation. In addition, this gene product induces mitochondria to release the apoptogenic proteins cytochrome c and caspase-9. Mutations in this gene cause combined oxidative phosphorylation deficiency 6 (COXPD6), a severe mitochondrial encephalomyopathy, as well as Cowchock syndrome, also known as X-linked recessive Charcot-Marie-Tooth disease-4 (CMTX-4), a disorder resulting in neuropathy, and axonal and motor-sensory defects with deafness and cognitive disability. Alternative splicing results in multiple transcript variants. A related pseudogene has been identified on chromosome 10. [provided by RefSeq, Aug 2015]

Locus ID:

9131