

Product datasheet for **SC213592**

FAM135A (NM_001105531) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	FAM135A (NM_001105531) Human 3' UTR Clone
Symbol:	FAM135A
Synonyms:	KIAA1411
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_001105531
Insert Size:	1283 bp



[View online »](#)

Insert Sequence:	>SC213592 3'UTR clone of NM_001105531 The sequence shown below is from the reference sequence of NM_001105531. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC CTGGTTGCTGCCCTCAAATATTTCCAATAGTATAAAAGCATTGTTAGCGACTGGACAATTACCTCATTC AACAAATGTTTCAAATAATGTATTATATTTAAATGTAGATGCTGATAAGTTCTAAGAAATATTTATACCT TTTTATATGGAAGATAATTTATATCATCCATGTTTAGTGCTTTTTAAACATCAACTTTACTTTCTAGGT AATGTGGCTGTGCAATATTTTTTAATTTTATCTTTTACTTTTCTATTACTTTTTCATATATTTTGCT ACCTAAGTATTTTCAAGTAACTTTAAGCCCATACCTGTGCTGATTGTTTATTATTGGCTTTCCACAAT TCTTACATCAGACTACATTATATTAGAGACCATTATTGCTAGAATAGCATGGGATTTAAATTTTCTAA TACTGGGGGTATTATTTAGTTAATTATAAATTTTCTTTTACATTTTACTGTGTTTTAACTGGAAATA AAATTATGGCTGTACAATATATTTTTGAAATCAACTTCTGTAGTTCTAAAATACAACTTATCATAAC AATCAAACCAGGTAGTTCATATAAAACAGTGTAATACAAGTTTTCTATAAAGTCATTACTGTTGCTTAA ACATATTTTCATGCCTATTTAAATATATTTTCTACTGGTGATTTCAACATTATTTCTCATACTGACTTTT ATTACTGGAAATGTTCTGTACATGTTGGCAGCAGATAAAGATTTTTGAATGTTTGAATGCCCTCTGCC TTGATTTGGTTGGAATTTTTGCTAAATTTGGTAATGTTGCTTGAACCTTTATGACTACATTTTCTTTAAC TTTTTTCATGGACTTCTTATATGTACATAATAATTAATGTTGAAATTTATGAAATACTTTTATGAAT TTAGATAATTTTTAAATATTGTTAAATTTTATTGAACTAAAAAGTAATGTAAATAAAATAATTCATGTT AAAGATGGAACAAAATAATTAACCTTTACATGTTTGGTGATACAGATGCAATGTTTTGATATATGGAG ATGTTGAGTCTTTTACTTAAAGGTGCTGAATAGCATTAAATCACTATTTTCTTTTCTGTTTT ACTTGTGAAAATAAAAAATGCACTAAGGTTGGGTAGAAGTTCTGTTGCACTCACTAATTGTGACAGACA GAGGTTTTTGTAAGTATTATTGTACAATTGATGCATGTTTATTTTACGTTGTTATTGCCTCTGGTG TTAATAAATGAACAAATGGCTATCTGGAGGAACAGCTACAA ACGCGTAAGCGGCCGCGGCATCTAGATTCGAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
Restriction Sites:	Sgfl-MluI
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_001105531.3</u>
Locus ID:	57579
MW:	51.1