

Product datasheet for **SC213599**

ZNF627 (NM_145295) Human 3' UTR Clone

Product data:

Product Type: 3' UTR Clones

Product Name: ZNF627 (NM_145295) Human 3' UTR Clone

Symbol: ZNF627

Mammalian Cell Selection: Neomycin

Vector: pMirTarget (PS100062)

ACCN: NM_145295

Insert Size: 1302 bp

Insert Sequence: >SC213599 3'UTR clone of NM_145295
The sequence shown below is from the reference sequence of NM_145295. The complete sequence of this clone may contain minor differences, such as SNPs.
Blue=Stop Codon **Red**=Cloning site

```
GGCAAGTTGGACGCCGCAAGATCCGCGAGATTCTCATTAAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAACGATCGCC
AACGCTTCAGTTGTCCCAGTTCTTTCATGAGCATGAAAGGAGTCACATAGAGAAACCCATGAAAGTAA
GAAATTTGGGAAAGCCTTCAGTCCTTTCTGTTTCTTTCAACTACGTGAAAGGATTCACAGTGGAGAAAG
ACCCTGTAAGATAATTGGCTTTAAATTACGAGAGACTTGTGATAGGACAGTAAACCTAGAGTTGGAGT
TGGATCTCTGGATTGTGTTATGTCAAGTGTGGTAGGTTAGGAAGTATTTCCAGAAATCCATCCATT
TGTGATTCATGATACAATTCACCAAGTAACCTATCTTACATGAGATTTCGGAAGTAAGTTAAGAAGGCAT
TAGTCATGGTTTGAAGCACCATACAGGGAGACAGCTGTGTGAATACAGGCTGTATGGACACTTGCTTC
CATCCCATTTTCTGCTTCTTTGGGTTGCCAATCAAGAGTATCCTCAAAACGACTTGACTTTAATTTTC
TCGGAGGTGATAGGCTTCCACACAGGTCTCCAGAAGCCCTGCATTGAATATCCATCCACACTTTGGTTT
TCCTTCAGACATTATTATGTCTGTACTAGGCAACTAATTCAGACTGTCCTGGTTGGGAATATTCTGTGA
TGCTCTGACTCCCTAGTCTGTAGACGGAATTGGCATAACGGTCTAATTTGTGTAGTAAGCACCTTTGTT
CATACTAGTAGTACTGTATTCTTGATTACGCTGATAGCTACCATGCTGCTGTCAAAACCAACCAGAG
GGGAGCTTGTCTTCTGCTGTAGTGTGAGTGGCCTCAGGAGCTTTGATGTGAGAAATGAGCA
CTTTCCTCTATCAGGAAATTTCAAGTGTTCCTGTTATTTCGTAGCTCAATGTAATGCCTCAGTTCATTT
TCAGTTGTTTGATTATATGTGACTAATATGTATTTTTATTCAAACAAGACTTCTGTACATGTTTCTTC
AAAACAGTTTATTAAGTGTCTTCAGTCTTGATTACATCAAGTTTATAATTTTGGCAAAATGTTAAGAC
ACTGTGAAGTCAGCGTTAACCATGTGCATACAACCTAAGGAATTTTCTCCTCATGTAATTTTACT
TTTCATGCTTATATAGTTTCACTTTTATCTTCATAGTAATTTCTCATCTACTCATAATACCAAAAGTT
AAGTCATGCTGTTTTGTGTGCTCTCTTGCTAAAGACCGCAGAGACCATACCTGTTGTCAAAGAGGGTGT
AATAAACTGTAATAATAATCATGACCACAAACCATAAAGGACTGTGTGATGTATCAGTAA
ACGCGTAAGCGGCCGCGGCATCTAGATTCTGAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
```



[View online »](#)

Restriction Sites:	Sgfl-Mlul
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	NM_145295.4
Summary:	May be involved in transcriptional regulation.[UniProtKB/Swiss-Prot Function]
Locus ID:	199692
MW:	49.6