

Product datasheet for **SC213606**

CMTM3 (NM_144601) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	CMTM3 (NM_144601) Human 3' UTR Clone
Symbol:	CMTM3
Synonyms:	BNAS2; CKLFSF3
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_144601
Insert Size:	1280 bp



[View online »](#)

Insert Sequence:	>SC213606 3'UTR clone of NM_144601 The sequence shown below is from the reference sequence of NM_144601. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC GAAGAGAATTCGACTCGGACTCTGACTGAAGGCCTGGCGGGTGCCTTGGCAACCTGAGCCACACAGGC CTCCACCCTGCGCTCAGAGGGTCTGCTGGCGTTGGAGCGGAGGCTTGGACTTCTGAGTTGCAGAGGG GGCTGCGGACACAGCAGGCCCTTACAGCCTCAGGTTCTGCTGAGCCAGCCTACCAGGCTTGGCCCT CAGCTCAGCACTGTTGACCACGCTGCGTATGAGGGCATCTTGGGTATCCCACTCCTTCTCCCCATTTCT GTCCACAGGCCCTTACGCCCTTTAACGTCTCTGCCAAAACCAGCACAAGGAGACAAAGCAGAGCCTTG TCTGTATCTGGGAGCAGGTGTTCCATGCTGCTAGGTGGCGGGGTCTTCTGTTTCACTAAC AGGAACAAAGACAGAAACCATGACAGGGCTGCCCGCCAGGCCCGGTGGGTTTGTCTGCACTTGGTGC TCCTGCCACACCAGCCACTTTGGTGACAATGACCCTTCCAAGAATCTTGGTTCAAGGAGCACCAGTT CCCTCTTCTATTCTGAAGCAGGGAGAAATTGACCTTTGCCTTGTGCGCCAGGAAGTGGGGCTCGGCACC CATAACTAACACCTCCCACTTGGAAACCATGTCTTCTGGGGGTGAGATGACCATTCTGGGTCTAAGA CTGTTTCAAAGAAGAGCTCATAGACTGACTGGTCCAGAAGACAGAGGGTACAACAGTGGCATCACAGTG ACAGTGTGATGGGAGCTGGGCGGGCCAGCCAAACCTCCTTCTTCTAGAGCCAGCCAGCAGGAG GAGTTCCTGGACCCTCAGGACAGTGAATTCAGACCTCAGGGCAGGTCTATGGGCCACTGCAGGAGAT GAGACCAGCCTTCTGTGTTACCTAACGATTTATACTGTGTATCTGTCTTTGATGGAATTTGTAACCT TTTATATTTTTTATGCAAAAGCAGCTTCTTAACAGATGGCATTTTCTGTGACTCTAGGCCTCACAAAA GAGCCAGAGTTCTGGACCCATGTTTGGAGCATTTGTAGCCTTATTCTTTCGCTGTGAATCTCTTACCC TGAAAAAAGCCATAATGAATTAAGCCAGACTGACCACTTGCTTGGAGTGTGTGCTTGAACCAACAGAG GCAATACTGTTGGGTATTGTATCAGGCTTCAGTACAACTGGTAACACCAATGTGGATCCTGACAGCTT TCAGTTTTAGCAAAAATACACGTGAAATCTGACTACCA ACGCGTAAGCGGCCGCGCATCTAGATTGAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA CGAGATTCGATTCCACCGCGCCTTCTATGAAAGG
Restriction Sites:	Sgfl-MluI
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_144601.5</u>
Summary:	This gene belongs to the chemokine-like factor gene superfamily, a novel family that is similar to the chemokine and the transmembrane 4 superfamilies of signaling molecules. This gene is one of several chemokine-like factor genes located in a cluster on chromosome 16. Alternatively spliced transcript variants containing different 5' UTRs, but encoding the same protein, have been identified. [provided by RefSeq, Jul 2008]
Locus ID:	123920
MW:	46.1