

Product datasheet for **SC213619**

SLC35F3 (NM_173508) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	SLC35F3 (NM_173508) Human 3' UTR Clone
Symbol:	SLC35F3
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_173508
Insert Size:	1298 bp
Insert Sequence:	>SC213619 3'UTR clone of NM_173508 The sequence shown below is from the reference sequence of NM_173508. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAACGATCGCC
AGAAGAGCCCGCCTTCCCTTCGCCCGCTAACACCACTCCTCTAGAACTCGGTGGTAATGACTGGGAGGT
CTATTCCTGCCGGGAGGAACCTCAGTTGGGTAAAGGTGTACATACCTGTACAGTTTTGGTCATCTGCGGT
AAGTTCATGGTATTTATTGGCATGTCCAATTTGCTTGCACTTTTCCACATCAGACCTTCCACTTAAGG
GTACCAATTCATACAAAAGAGAAAAGAAGAACCTCTCAGTGCTTTTCCAACGAATGTAAGTGGGTTT
AAAAGTTCCTGCGTAGATCGTTTTGGATGGCTGACGTTCTTGACAAGCCAGCCATCAGGGCAAGTGTT
TTGAATCTTAGCAAGTGTTTGAATCTTAGCAACTGAACATGACATTGCACACTGTGATTATTTTCAG
CTATGGTAGGTCAATTTTGTATATACCAGAGCTGTACTCTTAATTTTAAGGAATTTCAATGAACCAA
AAGATAACTGTCAATTTATTTGGATGGATGGATGGAGGGATGGAAGGAGGGAGGGAAGGATGTCTGGGA
TGTGAAGAGGCAAGGCAGAAATGGATGATGAACATCAAAAACCTGTCAGAGGACAGGTCTGTGCCAGGTA
CCAGGAAGGACAAGATAAGTCCAAGCCCCGGCTGGCCTGCTCTCACGTGTCTGTGAAGATCATGCAAG
AGAAATCACAGTGCCTTCTACAGAATTTTCGTCTTTACCTATGTGAAGCGAGGTGACGTGATACGTAC
TGGCGCCGCTTTATAATTTAGATGTAAAAATCTTTAGAAACAAATAAACTCTCTATATATGTGTATGT
CTGTGTACAAAAAATGACAGAGCTGATGGCCAGTGTATACAGAGCGTGGCCCGCGGTGTACAATACCC
ATATAAGGTACATTGTGCAGGAGGGGAATTGCTGGCTGCTTTTACTTCCTGACCAAGACTGAAAAATTA
TTTACTGAAATCTGTAAACCTTTTATGAACTTTTAAGCACCAGGCTGTTTACTTACACAATTTAGGT
CTGCCAGAAAAATCTATCTGTGATAGATCTGTAAAGAGGGTCAGGGGTTAGAGTTTACTATTTTGAAG
TTTACATTGTTACATATGAAATGGAAACATTATTTGAAACGTTGTACATAACCAATGGTGCATTCTGT
AACCATGGAGTCTTCTGTTTCTGGGGGAAAGGGGCATTCATGACCTGAACTTTTAGCAAATTATTAT
TCTCAGTTTCCATTACCTGTTTGGCCAAACAGATTAATAAAATATTTGAAAAAGAA
ACGCGTAAGCGGCCGCGCATCTAGATTGGAAGAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
  
```


[View online »](#)

Restriction Sites:	Sgfl-Mlul
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_173508.4</u>
Summary:	May be a thiamine transporter.[UniProtKB/Swiss-Prot Function]
Locus ID:	148641
MW:	50.3