

Product datasheet for **SC213641**

TMEM25 (NM_032780) Human 3' UTR Clone

Product data:

Product Type: 3' UTR Clones
Product Name: TMEM25 (NM_032780) Human 3' UTR Clone
Vector: pMirTarget (PS100062)
Symbol: TMEM25
ACCN: NM_032780
Insert Size: 1279 bp
Insert Sequence: >SC213641 3'UTR clone of NM_032780

The sequence shown below is from the reference sequence of NM_032780. The complete sequence of this clone may contain minor differences, such as SNPs.
Blue=Stop Codon **Red**=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAAACAATTGGCAGAGCTCAGAATTCAACGCATCGCC
AGCGTGAGCAGTGATGAGATCTGGCTCTGAGCCGAGGGCGAGACAGGAGTATTCTCTTGGCCTCTGGAC
ACCCTCCCGTTTCTCCAAGGCATCCTCTACCTAGCTAGGTACCAACGTGAAGAAGTTATGCCACTGCC
ACTTTTGCTTGGCCTCCTGGCTGGGGTGCCTCCATGTCATGCACGTGATGCATTTCACTGGGCTGTAA
CCCGCAGGGGCACAGGTATCTTTGGCAAGGCTACCAAGTTGGACGTAAGCCCTCATGCTGACTCAGGGT
GGGCCCTGCATGTGACTGGGCCCTTCCAGAGGGAGCTCTTTGGCCAGGGGTGTTAGATGTCATCC
AGCATCCAAGTGTGGCATGGCCTGCTGTATACCCACCCAGTACTCCACAGCACCTTGTACAGTAGGC
ATGGGGCGTGCCTGTGTGGGGACAGGGAGGGCCCTGCATGGATTTTCTCCTTCTATGCTATGTAG
CCTTGTTCCTCAGGTAATAATTTAGGACCCTGCTAGCTGTGCAGAACCAATTGCCCTTTCACAGAAA
CCAACCCTGACCCAGCGGTACCGGCCAAGCACAAACGTCCTTTTGTGTCACACGTCTCTGCCCTTCA
CTTCTTCTTCTGTCCCACCTCCTCTTGGGAATTCTAGGTTACACGTTGGACCTTCTCTACTACTTC
ACTGGGCACTAGACTTTTCTATTGGCCTGTGCCATCGCCAGTATTAGCACAAGTTAGGGAGGAAGAGG
CAGGCGATGAGTCTAGTAGCACCCAGGACGGCTTGTAGCTATGCATCATTTTCTACGGCGTTAGCACT
TTAAGCACATCCCTAGGGGAGGGGTGAGTGAGGGGCCAGAGCCCTCTTTGTGGCTTCCCCACGTTT
GGCCTTCTGGGATTCACTGTGAGTGTCTGAGCTCTCGGGGTTGATGGTTTTTCTCTCAGCATGTCTCC
TCCACCACGGGACCCAGCCCTGACCAACCCATGGTTGCCTCATCAGCAGGAAGGTGCCCTTCTGGAG
GATGGTCGCCACAGGCACATAATTCAACAGTGTGAAGCTTTAGGGGAACATGGAGAAAGAAGGAGACC
ACATACCCCAAAGTGACCTAAGAACAACCTTAAAAAGCAACATGTAATGATTGGAATTAATATAGTAC
AGAATATATTTTTCCCTTGTGAGATCTCTTTTGTAAATGTTTTTTCATGTTACTGCCTAGGGCGGTGCT
GAGCACACAGCAAGTTTAAATAAACTTGACTGAATTCA
ACGCGTAAGCGGCCGCGGCATCTAGATTGGAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
```

Restriction Sites: Sgfl-MluI



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	NM_032780.4
Summary:	In neurons, modulates the degradation of NMDA receptor GRIN2B subunit. Plays a role in the regulation of neuronal excitability.[UniProtKB/Swiss-Prot Function]
Locus ID:	84866
MW:	46.7