

Product datasheet for **SC213642**

TMEM25 (NM_001144034) Human 3' UTR Clone

Product data:

Product Type: 3' UTR Clones
Product Name: TMEM25 (NM_001144034) Human 3' UTR Clone
Vector: pMirTarget (PS100062)
Symbol: TMEM25
ACCN: NM_001144034
Insert Size: 1279 bp
Insert Sequence: >SC213642 3'UTR clone of NM_001144034

The sequence shown below is from the reference sequence of NM_001144034. The complete sequence of this clone may contain minor differences, such as SNPs.
 Blue=Stop Codon Red=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
AGCGTGAGCAGTGATGAGATCTGGCTCTGAGCCGAGGGCGAGACAGGAGTATTCTCTTGGCCTCTGGAC
ACCCTCCCGTTTCTCCAAGGCATCCTCTACCTAGCTAGGTACCAACGTGAAGAAGTTATGCCACTGCC
ACTTTTGCTTGGCCTCCTGGCTGGGGTGCCTCCATGTCATGCACGTGATGCATTTCACTGGGCTGTAA
CCCGCAGGGGCACAGGTATCTTTGGCAAGGCTACCAAGTTGGACGTAAGCCCTCATGCTGACTCAGGGT
GGGCCCTGCATGTGATGACTGGGCCCTCCAGAGGGAGCTCTTTGGCCAGGGGTGTTAGATGTCATCC
AGCATCCAAGTGTGGCATGGCCTGCTGTATACCCACCCAGTACTCCACAGCACCTTGTACAGTAGGC
ATGGGGCGTGCCTGTGTGGGGACAGGGAGGGCCCTGCATGGATTTTCTCCTTCTATGCTATGTAG
CCTTGTTCCTCAGGTAATAATTTAGGACCCTGCTAGCTGTGCAGAACCAATTGCCCTTTCACAGAAA
CCAACCCTGACCCAGCGGTACCGGCCAAGCACAAACGTCCTTTTGTGTCACACGTCTCTGCCCTTCA
CTTCTTCTTCTGTCCCACCTCCTCTTGGGAATTCTAGGTTACACGTTGGACCTTCTCTACTACTTC
ACTGGGCACTAGACTTTTCTATTGGCCTGTGCCATCGCCAGTATTAGCACAAGTTAGGGAGGAAGAGG
CAGGCGATGAGTCTAGTAGCACCCAGGACGGCTTGTAGCTATGCATCATTTCCTACGGCGTTAGCACT
TTAAGCACATCCCTAGGGGAGGGGTGAGTGAGGGGCCAGAGCCCTCTTTGTGGCTTCCCCACGTTT
GGCCTTCTGGGATTAAGTGTGAGTGTCTGAGCTCTCGGGTTGATGGTTTTTCTCTCAGCATGTCTCC
TCCACCACGGGACCCAGCCCTGACCAACCCATGGTTGCCTCATCAGCAGGAAGGTGCCCTTCTGGAG
GATGGTCGCCACAGGCACATAATTCAACAGTGTGAAGCTTTAGGGGAACATGGAGAAAGAAGGAGACC
ACATACCCCAAAGTGACCTAAGAACAACCTTAAAAAGCAACATGTAATGATTGGAATTAATATAGTAC
AGAATATATTTTTCCCTTGTGAGATCTCTTTTGTAAATGTTTTTTCATGTTACTGCCTAGGGCGGTGCT
GAGCACACAGCAAGTTTAAATAAACTTGACTGAATTCA
ACGCGTAAGCGGCCGCGCATCTAGATTGAAAGAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTCGATTCACCGCCGCTTCTATGAAAGG
```

Restriction Sites: SgfI-MluI



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	NM_001144034.2
Summary:	In neurons, modulates the degradation of NMDA receptor GRIN2B subunit. Plays a role in the regulation of neuronal excitability.[UniProtKB/Swiss-Prot Function]
Locus ID:	84866
MW:	46.7