

Product datasheet for **SC213702**

RABL2B (NM_007081) Human 3' UTR Clone

Product data:

Product Type: 3' UTR Clones

Product Name: RABL2B (NM_007081) Human 3' UTR Clone

Vector: pMirTarget (PS100062)

Symbol: RABL2B

ACCN: NM_007081

Insert Size: 1319 bp

Insert Sequence: >SC213702 3' UTR clone of NM_007081
The sequence shown below is from the reference sequence of NM_007081. The complete sequence of this clone may contain minor differences, such as SNPs. **Red**=Cloning site
Blue=Stop Codon

CAATTGGCAGAGCTCAGAATTCAAGCGATCGC

ATCGAGACCCCATCAGAGGAGGGCGCCTCTCCCCACAGCTGAGGGGCTGGGGCTAGGGGTGGGTGGAGCC
CTTTTAAATACCCCTTCCCTTCAACAACCTCTCCAGCTCTGAATGGAGAACTCTCTAGGCCATCCCTCT
TCTACCTCCTGCAACCCACCCATCCTATTAGCCTCCCACATTCTAGGCCGTGATACAGGGATGAGGTCA
GCACCAGCAAACTCTGGACTGGTGAAGAATCCCCACCAGATCTCCTTGAAGCAGAAGTAGGGATCAGC
ATCATTAACACCTTCCCCACCCCTCCCCGAGGCAGACAGTGAAGAGAATCAGAAAACATGATTATGTG
TCACTTAATACAGGAATTTAGGTGTTTTTGGTGTTTTTGTGTTTTGTTTTGTTTTTCTTTCCAAAGCT
CACCTCGGGGACAATTCCTTGGGCTTCTCCTGAGGTAATGATTTACCCCCACCCACAGCTGAGTCTGT
GAGGCCCATCCTTTCCCTACGTTTTCTCCATCTTTTTCTCTTTCAGTCTCCAGTCATCTGGTTTGT
TTGTTTCTTTGTTGTTCTGAGACGGAGTCTCGCTCTGTCGCCAGGCTGGAGTGCAGTGGCGCAGTCTCG
GCTCGCTGCAACCTCTGACTCCCTGGTTCAAACGATTCTCCTGCCTCAGCCTCCCGAGTGGCTGGCATCA
CCACGCCAGCTAATTTTTGTATTTTTAGTAGAGACGGGTTTACCATGTTGCCAGGATGGTCTCGAT
CTCCTCACCTCGTGATCCGCCCGCTCGGCCTCCCAAAGTGCTGGGATTACAGGCATGAGCCACCGCGCC
CGGCCCAATCATCTGTTTTTAAACAATCGTTTTTGAGCAGATAGCTATTCATTCCAGATTTTCGTGTAC
CCATTCTGTTTCAGGAGCTCTTCTAGGTAAAGCTGAGATCACAGGAACAGCAGGTGACAGGCCTAGCTAT
AGTTAGGAATACACAAGCCGTAATTCGAGTCTTACAGCCATACCACAAGGTACGTCCATTTGGACTAC
AAGAAGAGCTTCTTTAAAGTTCTATTTACAGATAAAGAGGCTGTCCTTTTTTTGAGGAATAGTTTGG
ACCTTGTCCTCCTGTGGGAGGCTGAGGACTGCAAGAGGAGAGCTAGCAGATATGCCTGTTCCACCCCTCT
CTGGTACTTGTGGCTTGTAGTATGTTTTATGATAATCTCGGCATTGTTGATTGTTTATTAATA
GGGTTTTGTTTTGTTGTTTCCTTTTTTACAGTAAAGGCTGAATTACATAAACATTGAG

ACGCGTAAGCGGCCGCGGCATCTAGATTCAAGAAAATGACCG

Restriction Sites: SgfI-MluI



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_007081.2</u>
Summary:	The RABL2B protein is a member of the RAB gene family which belongs to the RAS GTPase superfamily. RABL2B is located within a subtelomeric region of 22q13.3. Multiple alternatively spliced transcript variants encoding several different isoforms have been found for this gene. [provided by RefSeq, Aug 2008]
Locus ID:	11158