

## Product datasheet for **SC213704**

### **RABL2B (NM\_001130921) Human 3' UTR Clone**

#### Product data:

Product Type:	3' UTR Clones
Product Name:	RABL2B (NM_001130921) Human 3' UTR Clone
Vector:	pMirTarget (PS100062)
Symbol:	RABL2B
ACCN:	NM_001130921
Insert Size:	1319 bp
Insert Sequence:	>SC213704 3' UTR clone of NM_001130921 The sequence shown below is from the reference sequence of NM_001130921. The complete sequence of this clone may contain minor differences, such as SNPs. <b>Red</b> =Cloning site <b>Blue</b> =Stop Codon

CAATTGGCAGAGCTCAGAATTCAA**GCGATCGC**

ATCGAGACCCCATCAGAGGAGGGCGGCTCTCCACAGCT**G**AGGGGCTGGGGCTAGGGTGGGTGGAGCC  
CTTTAAATACCTTCCCTTCAACAACCTCCAGCTCTGAATGGAGAACTCTTAGGCCATCCCTCT  
TCTACCTCCTGCAACCCACCCATCCTATTAGCCTCCCACATTCTAGGCCGTGATACAGGGATGAGGTCA  
GCACCAGCAAACCTGGACTGGTGAAGAATCCCCACAGATCTCCTTGAAGCAGAAGTAGGGATCAGC  
ATCATTAAACACCTTCCCACCCCTCCCCGAGGCAGACAGTGAAGAGAATCAGAAAACATGATTATGTG  
TCACTTAATACAGGAAATTTAGGTGTTTTTGGTGTTTTTGTGTTTTGTGTTTTGTGTTTTCCAAAGCT  
CACCTCGGGGACAATTCCTTGGGCTTCTCCTGAGGTAATGATTACCCCCACCCACAGCTGAGTCTGT  
GAGGCCCATCCTTCCCTACGTTTTCTCCATCTTTTTCTCTTTCAGTCTCCAGTCATCTGGTTTTGT  
TTGTTTTCTTGTTCGTTCTGAGACGGAGTCTCGCTCTGTCGCCAGGCTGGAGTGCAGTGGCCAGTCTCG  
GCTCGCTGCAACCTCTGACTCCCTGGTTCAAACGATTCTCCTGCCTCAGCCTCCCGAGTGGCTGGCATCA  
CCACGCCAGCTAATTTTTGTATTTTTAGTAGAGACGGGTTTTACCATGTTGCCAGGATGGTCTCGAT  
CTCCTCACCTCGTATCCGCCCGCTCGGCTCCCAAAGTGTGGGATTACAGGCATGAGCCACCCGCC  
CGGCCCAATCATCTGTTTTAAACAATCGTTTTTGGAGCAGATAGCTATTCATTCCAGATTTTCGTGTAC  
CCATTCTGTTTCAGGAGCTCTTCTAGGTAAGCTGAGATCACAGGAACAGCAGGTGACAGGCCTAGCTAT  
AGTTAGGAATACACAAGCCGTAATAACGAGTCTTACAGCCATACCACAAGGTACGTCCATTTGGACTAC  
AAGAAGAGCTTCTTTAAAGTTCTATTTTACAGATAAAGAGGCTGTCTTTTTTTTGGAGGAATAGTTTGG  
ACTTGTGCCTCCTGTGGGAGGCTGAGGACTGCAAGAGGAGAGCTAGCAGATATGCCTGTTCCACCCCTCT  
CTGGTACTTGTGGCTGTAGTATGTTTTATGATAATCTCGGGCATTGTTTGCATTGTGTTTATAATA  
GGTTTTGTTTTGTTGTTTCTTTTTTACAGTAAAGGCTGAATTACATAAACATTGAG

**ACGCGT**AAGCGGCCGCGGCATCTAGATTCAAGAAAATGACCC

Restriction Sites: SgfI-MluI



[View online »](#)

<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
<b>Components:</b>	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
<b>RefSeq:</b>	<u><a href="#">NM_001130921.1</a></u>
<b>Summary:</b>	The RABL2B protein is a member of the RAB gene family which belongs to the RAS GTPase superfamily. RABL2B is located within a subtelomeric region of 22q13.3. Multiple alternatively spliced transcript variants encoding several different isoforms have been found for this gene. [provided by RefSeq, Aug 2008]
<b>Locus ID:</b>	11158