

## Product datasheet for **SC213754**

### ZNF135 (NM\_003436) Human 3' UTR Clone

#### Product data:

**Product Type:** 3' UTR Clones  
**Product Name:** ZNF135 (NM\_003436) Human 3' UTR Clone  
**Vector:** pMirTarget (PS100062)  
**Symbol:** ZNF135  
**Synonyms:** pHZ-17; pT3; ZNF61; ZNF78L1  
**ACCN:** NM\_003436  
**Insert Size:** 1324 bp  
**Insert Sequence:** >SC213754 3'UTR clone of NM\_003436

The sequence shown below is from the reference sequence of NM\_003436. The complete sequence of this clone may contain minor differences, such as SNPs.

Blue=Stop Codon Red=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
ACCAAGCACCAGAGAACTCACACTGGATTAACCCACTCCACATGTGCTGGGGACATAGGAAGACCTTAA
GCCATAGCTCATCCCTTTCTAGATTTGACCCAATCATAACACATGAGAAACGTACATTCATACACAAGCC
TTTTCACACAGCACTCCCTCAGACACCCCTCAGAGAGTTCACACTGATGGGAAATGACCATGGGACCAC
CAAGCTCTAGGTCATCCATCTCTGCATCCAAATAGTAGGGAAACGTGGAGATAATCAACTCAGGACC
TTCAGCCTTGAACGCCATTAGTGCTATGTTATAGAACCACAAAAAGAAATGGAACAAATGTAGTGG
ATCCAGGGAACGCTTTTGTCCAAGGATTCACCGTATTCCAACCAGAGATGTTCAAATTGGTGAGAAAC
CCAACAAATGCCTTTCATATACGAGAACCAATGAAGTCAGAAATTTGCCATTATTGCACATCACATT
TTTGGGGGGAAAGTGCTTATGAATGGTGCAGGTTGACTCTGATATTCATCCCAATGACAGTATGGC
AGAGTGTCCAGAAATGAGAGTGGCATCTTATGGAATCACTGTGGATACTGACTGTCTCAGTAAACAG
CTGTCTGTCTGTGTATATTGTTGATCAGGATACATGGCAGCCAGTCACAGATTGGAATTACATAT
GACAAAGTATCAGTGTACTATAAACAGTTTTTATGTTATCCCTGCATTATTTTGAATTAATCTTTAT
ATGCAATGAGATTGAAAAGCTTTGTATGGGAAGACTCAAATGTAAGCTGCTTCCGTAGAGTCTCACT
GATTCTGAGATGGCTTTTGTGCTCTGTTCTCCCTCTACATTTCTCTGCAGAACTCGCGTTAGAAACAC
AGATATTTGTTTTACAAAAGGGAGATTTTCTTTGTTAAACCATCGTCTTATAAGCAATAGCAAATT
CATGTTAGAACTTCTGTTTTGCAAGATTCATCTTTTTGGGATATGTTCAAGTTAATCCTCTCAAAC
GCTTTTTATTGTTTTACAAATTTAACATTTTGTAAACCCTAAGTCCACAAATAGCAGTCTTTCTG
ACAACATAACCTTTAAATGGTGACTTGGCTGCCCTCATTAGAAATGCATTGGCTGGCCAGGCGCGGTG
ACTTAGCCTGTAATCCCAGCACTTTGGGAGGCCAAGGTGGGAGGATCACCTGAGGTCGGGAGATGGAG
GTTGCAGTGAGCTGAGATCGTACCATTGCACTCCAGCCTGGGCAACAAGAGTGAACCTCTGTCTCAGAA
AAAAAAAAAAAAA
ACGCGTAAGCGGCCGCGCATCTAGATTCGAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCCTTCTATGAAAGG
```



[View online >](#)

<b>Restriction Sites:</b>	Sgfl-Mlul
<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
<b>Components:</b>	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
<b>RefSeq:</b>	<a href="#">NM_003436.4</a>
<b>Summary:</b>	Plays a role in the regulation of cell morphology and cytoskeletal organization. May be involved in transcriptional regulation.[UniProtKB/Swiss-Prot Function]
<b>Locus ID:</b>	7694
<b>MW:</b>	50.7