

Product datasheet for **SC213760**

GPR92 (LPAR5) (NM_001142961) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	GPR92 (LPAR5) (NM_001142961) Human 3' UTR Clone
Symbol:	GPR92
Synonyms:	GPR92; GPR93; KPG_010; LPA5
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_001142961
Insert Size:	1325 bp



[View online »](#)

Insert Sequence: >SC213760 3'UTR clone of NM_001142961
The sequence shown below is from the reference sequence of NM_001142961. The complete sequence of this clone may contain minor differences, such as SNPs.
Blue=Stop Codon **Red**=Cloning site

CGAAGTTGGACGCCCGCAAGATCCGCAGATTCTCATTAAAGGCCAAGAAGGCCGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAA**CGCATCGCC**
ACACAGTGTCCCCAGGATTCGCCCTCT**AG**ACACACATGCCATTGCGCTGTCCGTGCCCGACTCCCAAC
GCCTCTCGTTCTGGGAGGCTTACAGGGTGTACACACAAGAAGGTGGGCTGGGCACCTTGGACCTTTGGG
GGCAATTCAGCTTAGCAACGCAGAAGAGTACAAAGTGTGGAAGCCAGGGCCAGGGAAGGCAGTGCTG
CTGGAATAGCTTCTTTAACTGTGAGCACGCAGAGCACCCCTTCTCCAGCGGTGGGAAGTGATGCAGA
GAGCCCACCCGTGCAGAGGGCAGAAGAGGACGAAATGCCTTTGGGTGGCAGGGCATTAACTGCTAAA
AGCTGGTTAGATGGAACAGAAAATGGGCATTCTGGATCTAAACGCCACAGGGGCTGAGAGCTGAAGA
GCACCAGTTTGGTGGACAAAGCTACTGAGATGCCTGTTCATCTGCTGACTTCTGTCTAGGCTCATGGA
TGCCACCCCCTTTCATTTCCGCCTAGGCTTCCCCTGCTACCACTGAGGCCTAATACAAGAGTTCCTAT
GGACAGAACTACATTCTTCTCGCATAGTGACTTGTGACAATTTAGACTTGGCATCCAGCATGGGATAG
TTGGGGCAAGGCAAAACTAACTTAGAGTTTCCCCCTCAACAACATCCAAGTCCAAACCTTTTTAGGT
ATCCTTTCTTCCATCACATCCCCTTTCCAGGCCTCCTCCATTTTAGTCCCTAATATTCTTTCTTTT
CTCTCTCTCTCGTTTCTCTTCC
TCTCTCTCTTTGTCAGAGTAAGGATAAAATTCTTTCTACTAAAGCACTGGTCTCAAACCTTTTTGGT
CTCAGACCCCACTCTTAGAAATTGAGGATCTCAAAGAGCTTTGCTTATATTTGTTCTTTTGATACTTA
CCATACTAGAAATTAAGCGAATACATTTTTAAAAATAATACACATGCACACATTACATTAGCCATGGG
AGCAATAATGTCACCACACACACTTCATGAAGCCTCTGGAAAACCTCTACAGTATACTTGTGAGAGAATG
AGAGTGAAAGGGACAATAACATCTGTGTAGCAGTATTATGAAAATAGCTTGACCTCGTGGACTTCCTC
AGAGGGTTGGTCCCTGGATCACACTTTGAGAACCATACTTGTCTGAAGTATTGGAGTTCATGTCTAAC
TTCTTCCCAGGCATTATGTACAGTGCTTTTTATTACTGTGGGGAGAGGGCAGTGCTAAATAAATTAAT
CACTACTGATAGTC
ACCGGTAAGCGGCGCGGCATCTAGATTCGAAGAAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTTTGATTCCACCGCGCCTTCTATGAAAGG

Restriction Sites: SgfI-MluI

OTI Disclaimer: Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).

Components: The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.

RefSeq: NM 001142961.1

Summary: This gene encodes a member of the rhodopsin class of G protein-coupled transmembrane receptors. This protein transmits extracellular signals from lysophosphatidic acid to cells through heterotrimeric G proteins and mediates numerous cellular processes. Many G protein receptors serve as targets for pharmaceutical drugs. Transcript variants of this gene have been described.[provided by RefSeq, Dec 2008]

Locus ID: 57121

MW: 48.6